lnBody270

Great precision for flexible solution



Modern

Ergonomic design for efficient workflow

Reliable

Proven and reliable performance

Standard


The archetype of body composition test

The Evolution of Health Indicator

Paradigm shift of your health indicator

While simple and inexpensive, BMI often faced vigorous debate regarding health evaluation.

Medical professionals recognized that it is not enough to evaluate health by only considering weight and height. Therefore, people start to think inside of the body - Body Composition. Based on that, the concept of FFMI which means how much fat free mass you have inside of your body emerged and it is now considered as a new standard of health evaluation.



After using weight scale to measure human body, people started to think of weight and health together. BMI was the indicator of whether you are underweight, normal or overweigh by only using height and weight.

$$BMI = \frac{weight (kg)}{height^2 (m^2)}$$

PBF(Percent Body Fat)

Percent Body Fat (PBF) is the percentage of body fat compared to body weight. Actual obese status can be known through PBF.


FFMI (Fat Free Mass Index)

Fat Free Mass Index (FFMI) is calculated by dividing Fat Free Mass by height squared. This index can be used for monitoring FFM.

InBody Applications

InBody Test results can be used as an indicator for proper prescription by monitoring your health status in many areas such as obesity, nutrition, rehabilitation etc.

Providing a smart healthcare service with high technology and helping everyone live a healthier life, that's the great mission of InBody.











Timeless Style, Yet Significant Feature Changes

The all new InBody270 is now more elegant and premium than ever before. Its silhouette, combined with clean lines and high-class design elements provides a truly impressive appearance.

InBody270 Smart Usage Flow





^{*} Carrying Bag, Thermal Printer and Lookin' Body120 are optional items.



Barcode Reader


Quick member registration and ID number input available with barcode reader.

Foldable Structure

Foldable structure for easy installation and safe relocation.



Personal Data Input



Results Sheet

Abundant information printed by printer and thermal printer*.



User-friendly Interface ---

Easy operation through both touch screen and keypad.





Step-by-step posture guide for complete measurement.





InBody Results Sheet

Parameters on the right part of results sheet can be displayed as optional, depending on the customers' needs.

You can select and print Nutrition Evaluation, Body Balance Evaluation, Waist Circumference and Blood Pressure data in addition to items displayed in the results sheet.

1 Body Composition Analysis -----

Body weight is the sum of Total Body Water, Protein, Minerals, and Body Fat Mass.

Maintain a balanced body composition to stay healthy.

2 Muscle-Fat Analysis

Compare the bar lengths of Skeletal Muscle Mass and Body Fat Mass.

The longer the Skeletal Muscle Mass bar is compared to the Body Fat Mass bar, the stronger the body is.

3 Obesity Analysis

BMI is an index used to determine obesity by using height and weight. PBF is the percentage of body fat compared to body weight.

4 Segmental Lean Analysis

Shows how well developed your muscle is for each part of the body. See where you should work on more.

6 Segmental Fat Analysis

You can check and manage fat mass on each segmental part of the body. Monitor Percent Body Fat of each part and try to keep them in 'Normal' range.

6 Body Composition History -----

Track the history of the body compositional change. Take the InBody Test periodically to monitor your progress.

InBody

ID	Height	Age	Gender	Test
Jane Doe	156.9cm	51	Female	20

... Body Composition Analysis

Total amount of water in body	Total Body Water	(L)	26.5
For building muscles	Protein	(kg)	7.2
For strengthening bones	Minerals	(kg)	2.64
For storing excess energy	Body Fat Mass	(kg)	22.8
Sum of the above	Weight	(kg)	59.1

2 Muscle-Fat Analysis

		Uı	nder	N	lorma	1			Ov	er
Weight	(kg)	55	70	85	100	= 59	.1	145	160	
SMM Skeletal Muscle Mass	(kg)	70	80	90 1	9.6	110	120	130	140	1
Body Fat Mass	(kg)	40	60	80	100	160	220	22.8	340	4

3 Obesity Analysis

		U	nder		Norma	al l			٥١	/e
BMI Body Mass Index	(kg/m²)	10.0	15.0	18.5	21.0	^{25.0} 24	.0	35.0	40.0	4
PBF Percent Body Fat	(%)	8.0	13.0	18.0	23.0	28.0	33.0	38.0	43.0	4

4 Segmental Lean Analysis 5 Segmental Fat A

8		J	
1.71 kg		1.71 kg	1.7 kg
86.6 %		86.8 %	193.3 %
Normal	16.7 kg 93.6 %	Normal Z.	Over
Left	Normal	Right	Le
5.48 kg		5.40 kg	3.5 kg
87.9%		86.6 %	156.1 %
Under		Under	Normal

6 Body Composition History

Weight	65.3	63.9	62.4	61.8	62.3	60.9
SMM Skeletal Muscle Mass	20.1	20.0	19.7	19.7	19.8	19.7
PBF Percent Body Fat	41.3	40.7	39.2	39.0	39.4	38.6
▼ Recent □ Total	17.10.10 09:15	17.10.30 09:40	17.11.02 09:35	17.12.15 11:01	18.01.12 08:33	18.02.10 15:50

[InBody270]

Date & Time 8.04.03. 09 : 46

 $(26.4 \sim 32.2)$

 $(7.0 \sim 8.6)$

 $(2.44 \sim 2.98)$

 $(10.3 \sim 16.5)$

(43.9 ~ 59.5)

160 170

460 520

53.0 58.0

nalysis

Fat Mass

1.7 kg

192.0 % Over

3.4 kg 155.4 %

Normal

75 190 205

5.0 50.0 55.0

InBody

www.inbody.com

nBody Score

66/100 Points

* Total score that reflects the evaluation of body composition. A muscular person may score over

8 Weight Control

Target Weight 51.7 kg
Weight Control - 7.4 kg
Fat Control - 10.9 kg
Muscle Control + 3.5 kg

Obesity Evaluation

10 Waist-Hip Ratio

0.91

0.75 0.85

11 Visceral Fat Level

11

Low 10 High

PReasearch Parameters

Calorie Expenditure of Exercise

Golf	104	Gateball	112
Walking	118	Yoga	118
Badminton	134	Table Tennis	134
Tennis	177	Bicycling	177
Boxing	177	Racketball	177
Mountain Climbing	193	Jumping Rope	207
Aerobics	207	Jogging	207
Soccer	207	Swimming	207
Japanese Fencing	296	Racketball	296
Squash	296	Taekwondo	296
*Based on your curre			
*Based on 30 minute	duratio	n	

(B) Impedance

RA LA TR RL LL **Z**(Ω) 20 kHz 443.8 444.0 29.0 340.7 330.6 100 kHz 407.3 408.8 26.6 299.5 289.8

60.5 59.1 19.8 19.6 37.8 38.6

18.03.15 18.04.03 08:35 09:46

Copyright@1996~by InBody Co., Ltd. All rights reserved. BR-English-F3/230-B-131217

7 InBody Score

This score shows the evaluation of your body composition, which includes muscle, fat, and water in the body.

·· 8 Weight Control

See how your body measures up to the recommended Weight, Muscle Mass, and Body Fat Mass for a good balance. The '+' means to gain and the '-' means to lose.

Obesity Evaluation

Evaluates obesity based on BMI and Percent Body Fat.

Waist-Hip Ratio (WHR)

Waist-Hip Ratio is the ratio of waist circumference to hip circumference.

1 Visceral Fat Level

Visceral Fat Level is an indicator based on the estimated amount of fat surrounding internal organs in the abdomen.

12 Research Parameters


Various nutritional outputs are provided such as Fat Free Mass, Basal Metabolic Rate, Obesity Degree, Recommended calorie intake per day and more. To see a complete list, please scan the results interpretation QR code.

· 🔞 Impedance

Impedance is the resistance value measured when electrical currents are applied throughout the body. Based on the measured data, key body composition outputs can be analyzed. Impedance is also used for many research purposes.

The InBody Results Sheet for Children

Specially designed results sheet with Growth Graph is available for Children.



Thermal Results Sheet

You can select and print Total Body Water, Protein, Mineral and Waist Circumference data in addition to items displayed in the thermal results sheet.

InBody 04/03/2018 09:46
ID : Jane Doe Height : 156.9cm Age : 51 Gender : Female Weight : 59,1kg
Weight 59.1 kg
Muscle Mass 19.6 kg Mass of muscle attached to your bones. (Skeletal Muscle Mass)
Body Fat Mass 22.8 kg
Percent Body Fat 38.6 % Reference Range : Male adult 10~20% Female adult 18~28%
BMI 24.0 kg/m² Reference Range: Adult 18.5~25.0 kg/m²
Basal Metabolic Rate 1154 kcal Minimum number of calories needed to sustain life at a resting state.
Waist Hip Ratio 0.91 Reference Range : Male adult 0.80~0.90 Female adult 0.75~0.85
Visceral Fat Level 11 Reference Range :1~9
Segmental Lean Analysis
Under Under Segmental Fat Analysis
1.7 kg 1.7 kg 192.0 % Over 11.9 kg Over 239.8 % Over 를
3.5 kg 156.1 % Normal * Segmental fat is estimated.
InBody Score 66 point
Fat Control - 10.9 kg
Muscle Control + 3.5 kg
443.8 444.0 29.0 340.7 330.6 407.3 408.8 26.6 299.5 289.8
InBody TEL:02-501-3939 FAX:02-501-3978


InBody Technology

The InBody has revolutionized the field of BIA (Bioelectrical Impedance Analysis) and is currently the most advanced product in the market. With a 98% correlation with gold standard body composition methods such as DEXA, the InBody is trusted by top hospitals, gyms, and professional sports teams because of its precision. The InBody is fast, accurate, and non-invasive.



DIRECT SEGMENTAL MEASUREMENT

Because the trunk contains the majority of lean body mass, an independent trunk measurement is required for an accurate measurement of muscle mass. InBody measures the impedance for each segment of the body directly.



SIMULTANEOUS MULTI-FREQUENCY (SMF-BIA)

Utilizing InBody's latest advancement in BIA, the patented SMF-BIA, InBody accurately measures ICW and ECW at the same time, reducing variances in your results from the water flow for a faster snapshot of your body composition.



8-POINT TACTILE ELECTRODE

8-point tactile electrodes with thumb technology contributes to the high reproducibility of InBody results.



NO USE OF EMPIRICAL ESTIMATIONS

The InBody only uses impedance directly acquired from each subject allowing the InBody to always produce correct results without use of empirical estimations, such as gender and age.

Cross-Reference









	InBodyS10	InBody770	InBody570	InBody370
Frequency	1, 5, 50, 250, 500, 1000kHz	1, 5, 50, 250, 500, 1000kHz	5, 50, 500kHz	5, 50, 250kHz
Results Sheet Ouputs				'
Total Body Water	•	•	•	•
Protein	•	•	•	•
Soft Lean Mass	•	•	•	•
Minerals	•	•	•	•
Fat Free Mass	•	•	•	•
Body Fat Mass	•	•	•	•
Weight	•	•	•	•
Skeletal Muscle Mass	•	•	•	•
Body Mass Index	•	•	•	•
Percent Body Fat	•	•	•	•
Segmental Lean Analysis	•	•	•	•
Segmental Fat Analysis		•	•	•
Segmental ECW Ratio	•	•		
ECW Ratio Analysis	•	•	•	
Body Composition History	•	•	•	
InBody Score		•	•	•
Body Type		•	•	
Weight Control		•	•	•
Nutrition Evaluation		•	•	•
Obesity Evaluation		•	•	•
Body Balance Evaluation		•	•	•
Segmental Body Water Analysis	•	•		
Segmental ICW Analysis		•		
Segmental ECW Analysis		•		
Segmental Circumference		•	•	
Intracellular Water	•	•	•	
Extracellular Water	•	•	•	
Basal Metabolic Rate	•	•	•	•
Waist-Hip Ratio		•	•	•
Waist Circumference	•	•	•	
Visceral Fat Level		•	•	
Visceral Fat Area	•	•		
Obesity Degree		•	•	
Bone Mineral Content	•	•	•	
Body Cell Mass	•	•	•	
Arm Circumference	•	•	•	
Arm Muscle Circumference	•	•	•	
TBW/FFM	•	•	- -	
Fat Free Mass Index		•		
Fat Mass Index		•		
Recommended Calorie Intake				
Calorie Expenditure of Exercise				•
QR Code		•	•	
Reactance	•	•	- -	
Phase Angle		•		
Impedance		•	•	•
impodunce	•	Ū		

) *Bajos

InBody270

20, 100kHz
·
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•
•

Key Specifications

Bioelectrical	Bioelectrical 10 Impedance Measurements by Using 2 Differen				
Impedance Analysis (BIA)	Impedance (Z) Frequencies at Each of 5 Segments (Right Arm,				
Measurement Items		Left Arm, Trunk, Right Leg, and Left Leg)			
Electrode Method	Tetrapolar 8-Point Tactile Electrodes				
Measurement Method	Direct Segmental Multi-frequency Bioelectrical Impedance Analysis Method (DSM-BIA), Simultaneous Multi-frequency Impedance Measurement (SMF-BIA)				
Body Composition Calculation Method	No Empirical Estimation				

Feature Specifications

Optional Equipment	Stadiometer from InBody, Blood pressure monitor from InBody,
Optional Equipment	Thermal printer, SD400
I	Name, Address and Contact information can be shown on the
Logo	InBody Results Sheet.
Displaying Method	LCD Monitor, Data management software Lookin' Body
Tropa of maguita about	InBody Test results sheet, InBody Test results sheet for children,
Type of results sheet	Thermal results sheet
William C. Haman	Audible indication for test in progress, test complete, and
Voice Guidance	successfully saved settings changes
D. (.1	Test results can be saved if the member ID is utilized. Up to
Database	100,000 results can be saved.
Test Mode	Self Mode, Professional Mode
A durinintunt on Manne	Setup: Configure settings and manage data
Administrator Menu	Troubleshooting: Additional information to help use the InBody
LICD Thumb duine	Copy, backup, or restore the InBody test data (data can be viewed
USB Thumb drive	on Excel or Lookin' Body data management software)
Danca da Dandan	The member ID will be automatically inputted when the barcode
Barcode Reader	ID is scanned.
Daglaum data	Backup data saved in the InBody by using a USB Thumb Drive,
Backup data	Restore results on the InBody from a backup file.

Other Specifications

Applied Rating Current	200μΑ (±40μΑ)	
Adapter	Manufacturer	BridgePower Corp.
	Model	BPM040S12F07
	Power Input	AC 100 ~240V, 50/60Hz, 1.2A
	Power Output	DC 12V, 3.4A
Display Type	600 x 1024 7inch Color TFT LCD	
Internal Interface	Touchscreen, Keypad	
External Interface	RS-232C 1EA, USB HOST 2EA, USB SLAVE 1EA, LAN (10T)	
	1EA, Bluetooth 1EA, Wi-Fi 1EA	
Compatible Printer	Laser/Inkjet Printers recommended by InBody	
	* A list of printers compatible with the InBody can be found at	
	http://www.inbodyservice.com	
Dimension	356 (W) × 796 (L) × 995 (H): mm	
	14.0 (W) × 31.3 (L) × 39.2 (H) : inch	
Equipment Weight	14kg (30.9lbs)	
Testing Time	About 15 seconds	
Operation Environment	$10 \sim 40$ °C (50 ~ 104°F), 30 ~ 75% RH, 70 ~ 106kPa	
Storage Environment	$-10 \sim 70 ^{\circ}\text{C} (14 \sim 158 ^{\circ}\text{F}), 10 \sim 80\% \text{RH}, 50 \sim 106 ^{\circ}\text{kPa} (\text{No})$	
	Condensation)	
Testing Weight Range	10 ~ 250kg (22.0 ~ 551.21bs)	
Testing Age Range	3 ~ 99 years	
Height Range	95 ~ 220cm (3ft. 1.4in. ~ 7ft. 2.61in.)	

^{*} Specifications may change without prior notice.

InBody

InBody Co., Ltd. [Head Office]

TEL: +82-2-501-3939 FAX: +82-2-578-2716 Website: http://www.inbody.com E-mail: info@inbody.com InBody USA. [USA]

TEL: +1-323-932-6503 FAX: +1-323-952-5009 Website: http://www.inbodyusa.com E-mail: info@inbodyusa.com InBody Japan Inc. [JAPAN]

TEL: +81-03-5298-7667 FAX: +81-03-5298-7668 Website: http://www.inbody.co.jp E-mail: inbody@inbody.co.jp InBody China. [CHINA]

TEL: +86-21-64439738, 9739, 9705 FAX: +86-21-64439706

Website: http://www.inbodychina.com E-mail: sales@inbodychina.com

InBody is a total healthcare device manufacturer that has acquired over 80 patent rights across the globe.















