*...going one step further*

**G70**  
(1008494)

Latin	3
English	4
Deutsch	6
Español	8
Français	10
Português	12
Italiano	14
日本語	16
Русский	18
英文	20
Türkçe	22

**Latin**

- 1 Pulmo dexter
- 2 Pulmo dexter, lobus superior
- 3 Pulmo dexter, lobus superior, segmentum apicale [S I]
- 4 Pulmo dexter, lobus superior, bronchus segmentalis apicalis [B I]
- 5 Pulmo dexter, lobus superior, segmentum posterius [S II]
- 6 Pulmo dexter, lobus superior, bronchus segmentalis posterior [B II]
- 7 Pulmo dexter, lobus superior, segmentum anterius [S III]
- 8 Pulmo dexter, lobus superior, bronchus segmentalis anterior [B III]
- 9 Fissura horizontalis pulmonis dextri
- 10 Pulmo dexter, lobus medius
- 11 Pulmo dexter, lobus medius, segmentum laterale [S IV]
- 12 Pulmo dexter, lobus medius, bronchus segmentalis lateralis [B IV]
- 13 Pulmo dexter, lobus medius, segmentum mediale [S V]
- 14 Pulmo dexter, lobus medius, bronchus segmentalis medialis [B V]
- 15 Fissura obliqua pulmonis dextri
- 16 Pulmo dexter, lobus inferior
- 17 Pulmo dexter, lobus inferior, segmentum superius [S VI]
- 18 Pulmo dexter, lobus inferior, bronchus segmentalis superior [B VI]
- 19 Pulmo dexter, lobus inferior, segmentum basale mediale [segmentum cardiacum] [S VII]
- 20 Pulmo dexter, lobus inferior, bronchus segmentalis basalis medialis [bronchus cardiacus] [B VII]
- 21 Pulmo dexter, lobus inferior, segmentum basale anterius [S VIII]
- 22 Pulmo dexter, lobus inferior, bronchus segmentalis basalis anterior [B VIII]
- 23 Pulmo dexter, lobus inferior, segmentum basale laterale [S IX]
- 24 Pulmo dexter, lobus inferior, bronchus segmentalis basalis lateralis [B IX]
- 25 Pulmo dexter, lobus inferior, segmentum basale posterius [S X]
- 26 Pulmo dexter, lobus inferior, bronchus segmentalis basalis posterior [B X]
- 27 Pulmo sinister
- 28 Pulmo sinister, lobus superior
- 29 Pulmo sinister, lobus superior, segmentum apicoposterius [S I+II]
- 30 Pulmo sinister, lobus superior, bronchus segmentalis apicoposterius [B I+II]

**Latin**

- 31 Pulmo sinister, lobus superior, segmentum anterius [S III]
- 32 Pulmo sinister, lobus superior, bronchus segmentalis anterior [B III]
- 33 Pulmo sinister, lobus superior, segmentum lingulare superius [S IV]
- 34 Pulmo sinister, lobus superior, bronchus lingularis superior [B IV]
- 35 Pulmo sinister, lobus superior, segmentum lingulare inferius [S V]
- 36 Pulmo sinister, lobus superior, bronchus lingularis inferior [B V]
- 37 Fissura obliqua pulmonis sinistri
- 38 Pulmo sinister, lobus inferior
- 39 Pulmo sinister, lobus inferior, segmentum superius [S VI]
- 40 Pulmo sinister, lobus inferior, bronchus segmentalis superior [B VI]
- 41 Pulmo sinister, lobus inferior, segmentum basale anteromediale [S VII+VIII]
- 42 Pulmo sinister, lobus inferior, bronchus segmentalis basalis anteromedialis [B VII+VIII]
- 43 Pulmo sinister, lobus inferior, segmentum basale laterale [S IX]
- 44 Pulmo sinister, lobus inferior, bronchus segmentalis basalis lateralis [B IX]
- 45 Pulmo sinister, lobus inferior, segmentum basale posterius [S X]
- 46 Pulmo sinister, lobus inferior, bronchus segmentalis basalis posterior [B X]
- 47 Oesophagus
- 48 Trachea
- 49 Truncus brachiocephalicus
- 50 V. brachiocephalica sinistra
- 51 V. brachiocephalica dextra
- 52 A. pulmonalis dextra
- 53 Vv. pulmonales dextrae
- 54 V. cava superior
- 55 Aorta, pars ascendens
- 56 Cor, atrium dextrum
- 57 Cor, auricula dextra
- 58 A. coronaria dextra
- 59 Cor, ventriculus dexter
- 60 Pleura parietalis
- 61 Diaphragma
- 62 A. subclavia sinistra
- 63 A. carotis communis sinistra
- 64 Aorta, arcus [Arcus aortae]
- 65 Bronchus principalis sinister
- 66 A. pulmonalis sinistra
- 67 V. pulmonalis sinistra superior
- 68 Truncus pulmonalis
- 69 Bronchus principalis dexter
- 70 Cor, auricula sinistra
- 71 A. coronaria sinistra
- 72 Cor, ventriculus sinister
- 73 Cor, apex [Apex cordis]
- 74 V. cava inferior
- 75 Pericardium
- 76 Aorta, pars descendens
- 77 Medulla spinalis
- 78 Costa
- 79 M. latissimus dorsi
- 80 Vertebra thoracica
- 81 M. trapezius
- 82 M. erector spinae
- 83 V. azygos

# Segmented Lung, 20 parts

English

This “Thoracic Organs with Segmented Lungs” model shows the organs and structures of the chest cavity. These include the right and left lungs and the pulmonary lobes, heart, pericardium, the large vessels around the heart and the trachea with oesophagus. The lung segments can be taken apart; inside the lungs, the segmental bronchi are shown in different colours. Located under the chest cavity structure are the diaphragm and the associated section of the pleura. The anatomical structures behind the oesophagus and the thoracic aorta have been cut horizontally at the level of the lower thoracic spine.

As well as depicting the organs and structures of the thoracic cavity, this model helps to provide a better understanding of the lung segments as isolated components of the individual lobes, each with its own segmental bronchus. It is very important to be able to examine and identify all segments of the anatomy using bronchoscope and radiological methods as any section may become diseased requiring diagnosis and specifically targeted treatment or surgery.

- 1 Right lung
- 2 Right lung, superior lobe
- 3 Right lung, superior lobe, apical segment [S I]
- 4 Right lung, superior lobe, apical segmental bronchus [B I]
- 5 Right lung, superior lobe, posterior segment [S II]
- 6 Right lung, superior lobe, posterior segmental bronchus [B II]
- 7 Right lung, superior lobe, anterior segment [S III]
- 8 Right lung, superior lobe, anterior segmental bronchus [B III]
- 9 Horizontal fissure of right lung
- 10 Right lung, middle lobe
- 11 Right lung, middle lobe, lateral segment [S IV]
- 12 Right lung, middle lobe, lateral segmental bronchus [B IV]
- 13 Right lung, middle lobe, medial segment [S V]
- 14 Right lung, middle lobe, medial segmental bronchus [B V]
- 15 Oblique fissure of right lung
- 16 Right lung, inferior lobe
- 17 Right lung, inferior lobe, superior segment [S VI]
- 18 Right lung, inferior lobe, superior segmental bronchus [B VI]
- 19 Right lung, inferior lobe, medial basal segment [S VII]
- 20 Right lung, inferior lobe, medial basal segmental bronchus [B VII]
- 21 Right lung, inferior lobe, anterior basal segment [S VIII]
- 22 Right lung, inferior lobe, anterior basal segmental bronchus [B VIII]
- 23 Right lung, inferior lobe, lateral basal segment [S IX]
- 24 Right lung, inferior lobe, lateral basal segmental bronchus [B IX]
- 25 Right lung, inferior lobe, posterior basal segment [S X]
- 26 Right lung, inferior lobe, posterior basal segmental bronchus [B X]
- 27 Left lung
- 28 Left lung, superior lobe
- 29 Left lung, superior lobe, apicoposterior segment [S I+II]
- 30 Left lung, superior lobe, apicoposterior segmental bronchus [B I+II]
- 31 Left lung, superior lobe, anterior segment [S III]
- 32 Left lung, superior lobe, anterior segmental bronchus [B III]
- 33 Left lung, superior lobe, superior lingular segment [S IV]
- 34 Left lung, superior lobe, superior lingular bronchus [B IV]
- 35 Left lung, superior lobe, inferior lingular segment [S V]
- 36 Left lung, superior lobe, inferior lingular bronchus [B V]
- 37 Oblique fissure of left lung
- 38 Left lung, inferior lobe
- 39 Left lung, inferior lobe, superior segment [S VI]
- 40 Left lung, inferior lobe, superior segmental bronchus [B VI]
- 41 Left lung, inferior lobe, anteromedial basal segment [S VII+VIII]
- 42 Left lung, inferior lobe, anteromedial basal segmental bronchus [B VII+VIII]
- 43 Left lung, inferior lobe, lateral basal segment [S IX]
- 44 Left lung, inferior lobe, lateral basal segmental bronchus [B IX]
- 45 Left lung, inferior lobe, posterior basal segment [S X]
- 46 Left lung, inferior lobe, posterior basal segmental bronchus [B X]

English

- 47 Oesophagus
- 48 Trachea
- 49 Brachiocephalic trunk
- 50 Left brachiocephalic vein
- 51 Right brachiocephalic vein
- 52 Right pulmonary artery
- 53 Right pulmonary veins
- 54 Superior vena cava
- 55 Ascending aorta
- 56 Heart, right atrium
- 57 Heart, right auricle
- 58 Right coronary artery
- 59 Heart, right ventricle
- 60 Parietal pleura
- 61 Diaphragm
- 62 Left subclavian artery
- 63 Left common carotid artery
- 64 Arch of aorta; Aortic arch
- 65 Left main bronchus
- 66 Left pulmonary artery
- 67 Left pulmonary vein
- 68 Pulmonary trunk
- 69 Right main bronchus
- 70 Heart, left auricle
- 71 Left coronary artery
- 72 Heart, left ventricle
- 73 Apex of heart
- 74 Inferior vena cava
- 75 Pericardium
- 76 Descending aorta
- 77 Spinal cord
- 78 Rib
- 79 Latissimus dorsi
- 80 Thoracic vertebra
- 81 Trapezius
- 82 Erector spinae
- 83 Azygos vein

# Assembly Instructions

English

## Heart

1. Insert the heart as shown in the diagram on page 24.

## Right Lung (Green)

Position the model and arrange the segments as shown in the diagram on page 25.

2. First fit the medial basal segment [S VII] (19) of the inferior lobe onto the base on the inside, bring together the colour-edged surfaces as shown in the diagram on page 25. The correct position can be seen on the diagram on page 26. Then rotate the model anticlockwise so that it is in the position shown on the diagram on page 27.
3. Fit the posterior basal segment [section X] (25) of the inferior lobe onto the base at the back. Bring together the colour-edged surfaces shown on the diagram on page 27. The correct position can be seen on the diagram on page 28. Then rotate the model clockwise so that it is in roughly the position shown on the diagram on page 29.
4. Fit the lateral basal segment [section IX] (23) of the inferior lobe onto the base on the outside. Bring together the colour-edged surfaces shown on the diagram on page 29. The correct position can be seen on the diagram on page 30.
5. Now fit the anterior basal segment [section VIII] (21) of the inferior lobe onto the base at the front. Bring together the colour-edged surfaces shown on the diagram on page 30. The correct position can be seen on the diagram on page 31. Then rotate the model clockwise so that it is in roughly the position shown on the diagram on page 32.
6. Now fit the medial segment [S V] (13) of the middle lobe onto the inside. Bring together the colour-edged surfaces shown on the diagram on page 32. The correct position can be seen on the diagram on page 33. Then rotate the model anticlockwise so that it is in roughly the position shown on the diagram on page 34.
7. Now fit the lateral segment [S IV] (11) of the middle lobe onto the outside. Bring together the colour-edged surfaces shown on the diagram on page 34. The correct position can be seen on the diagram on page 35. Then rotate the model anticlockwise so that it is in roughly the position shown on the diagram on page 36.
8. Now fit the superior segment [S VI] (17) of the inferior lobe onto the top. Bring together the colour-edged surfaces shown on the diagram on page 36. The correct position can be seen on the diagram on page 37. Then rotate the model clockwise so that it is in roughly the position shown on the diagram on page 37.
9. Now fit the posterior segment [S II] (5) of the superior lobe onto the back. Bring together the colour-edged surfaces shown on the diagram on page 37. The correct position can be seen on the diagram on page 38. Then rotate the model clockwise so that it is in roughly the position shown on the diagram on page 38.
10. Now fit the anterior segment [S III] (7) of the superior lobe onto the front. Bring together the colour-edged surfaces shown on the diagram on page 38. The correct position can be seen on the left diagram on page 39.
11. Now fit the apical segment [S I] (3) of the superior lobe onto the top. Bring together the colour-edged surfaces shown on the left diagram on page 39. The correct position can be seen on the right diagram on page 39.

English

## Left Lung (Red)

Position the model and arrange the segments as shown in the diagram on page 40.

12. First fit the anteromedial basal segment [S VII + VIII] (41) of the inferior lobe onto the base at the front on the inside. Bring together the colour-edged surfaces shown on the diagram on page 40. The correct position can be seen on the diagram on page 41.
13. Then fit the lateral basal segment [section IX] (43) of the inferior lobe onto the base on the outside. Bring together the colour-edged surfaces shown on the diagram on page 42. The correct position can be seen on the diagram on page 43. Then rotate the model clockwise so that it is in roughly the position shown on the diagram on page 43.
14. Then fit the posterior basal segment [section X] (45) of the inferior lobe onto the base at the back. Bring together the colour-edged surfaces shown on the diagram on page 43. The correct position can be seen on the diagram on page 44.
15. Now fit the superior segment [S VI] (39) of the inferior lobe onto the top. Bring together the colour-edged surfaces shown on the diagram on page 45. The correct position can be seen on the diagram on page 46. Then rotate the model anticlockwise so that it is in roughly the position shown on the diagram on page 46.
16. Now fit the inferior lingular segment [S V] (35) of the superior lobe onto the bottom. Bring together the colour-edged surfaces shown on the diagram on page 46. The correct position can be seen on the left diagram on page 47.
17. Now fit the superior lingular segment [S IV] (33) of the superior lobe onto the top. Bring together the colour-edged surfaces shown on the left diagram on page 47. The correct position can be seen on the right diagram on page 47.
18. Now fit the anterior segment [S III] (31) of the superior lobe onto the front. Bring together the colour-edged surfaces shown on the right diagram on page 47. The correct position can be seen on the diagram on page 48. Then rotate the model clockwise so that it is roughly in the position shown on the left diagram on page 48.
19. Now fit the apicoposterior segment [S I + II] (29) of the superior lobe onto the top at the back. Bring together the colour-edged surfaces shown on the left diagram on page 48. The correct position can be seen on the right diagram on page 48.

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- |  |   |
|--|---|
| 1 Rechte Lunge   | 24 Rechte Lunge, Unterlappen, Außenbasissegmentbronchus [B IX]          |
| 2 Rechte Lunge, Oberlappen   | 25 Rechte Lunge, Unterlappen, Hinterbasissegment [S X]                  |
| 3 Rechte Lunge, Oberlappen, Spitzensegment [S I]                               | 26 Rechte Lunge, Unterlappen, Hinterbasissegmentbronchus [B X]          |
| 4 Rechte Lunge, Oberlappen, Spitzensegmentbronchus [B I]                       | 27 Linke Lunge  |
| 5 Rechte Lunge, Oberlappen, Hintersegment [S II]                               | 28 Linke Lunge, Oberlappen  |
| 6 Rechte Lunge, Oberlappen, Hintersegmentbronchus [B II]                       | 29 Linke Lunge, Oberlappen, Spitzen- und Hintersegment [S I+II]         |
| 7 Rechte Lunge, Oberlappen, Vordersegment [S III]                              | 30 Linke Lunge, Oberlappen, Spitzen- und Hintersegmentbronchus [B I+II] |
| 8 Rechte Lunge, Oberlappen, Vordersegmentbronchus [B III]                      | 31 Linke Lunge, Oberlappen, Vordersegment [S III]                       |
| 9 Horizontale Furche der rechten Lunge   | 32 Linke Lunge, Oberlappen, Vordersegmentbronchus [B III]               |
| 10 Rechte Lunge, Mittellappen  | 33 Linke Lunge, Oberlappen, Oberzungensegment [S IV]                    |
| 11 Rechte Lunge, Mittellappen, Außensegment [S IV]                             | 34 Linke Lunge, Oberlappen, Oberzungensegmentbronchus [B IV]            |
| 12 Rechte Lunge, Mittellappen, Außensegmentbronchus [B IV]                     | 35 Linke Lunge, Oberlappen, Unterzungensegment [S V]                    |
| 13 Rechte Lunge, Mittellappen, Innensegment [S V]                              | 36 Linke Lunge, Oberlappen, Unterzungensegmentbronchus [B V]            |
| 14 Rechte Lunge, Mittellappen, Innensegmentbronchus [B V]                      | 37 Schräge Furche der linken Lunge                                      |
| 15 Schräge Furche der rechten Lunge  | 38 Linke Lunge, Unterlappen   |
| 16 Rechte Lunge, Unterlappen   | 39 Linke Lunge, Unterlappen, Obersegment [S VI]                         |
| 17 Rechte Lunge, Unterlappen, Obersegment [S VI]                               | 40 Linke Lunge, Unterlappen, Obersegmentbronchus [B VI]                 |
| 18 Rechte Lunge, Unterlappen, Obersegmentbronchus [B VI]                       | 41 Linke Lunge, Unterlappen, Innen- und Vorderbasissegment [S VII+VIII] |
| 19 Rechte Lunge, Unterlappen, Innenbasissegment [S VII]                        | 42 Innen- und Vorderbasissegmentbronchus [B VII+VIII]                   |
| 20 Rechte Lunge, Unterlappen, Innenbasissegmentbronchus [Herzbronchus] [B VII] | 43 Linke Lunge, Unterlappen, Außenbasissegment [S IX]                   |
| 21 Rechte Lunge, Unterlappen, Vorderbasissegment [S VIII]                      | 44 Linke Lunge, Unterlappen, Außenbasissegmentbronchus [B IX]           |
| 22 Rechte Lunge, Unterlappen, Vorderbasissegmentbronchus [B VIII]              | 45 Linke Lunge, Unterlappen, Hinterbasissegment [S X]                   |
| 23 Rechte Lunge, Unterlappen, Außenbasissegment [S IX]                         |   |

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© Copyright 2013 for instruction manual and design of product:
3B Scientific GmbH
3B Scientific GmbH, Germany