ARTICOLI & APPARECCHI PER MEDICINA

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Our hot air sterilizers are designed to allow a quick and safe sterilization of surgical instruments without envelope, in accordance with denominated method dry heat. The sterilization chamber is entirely made of stainless steel resistant to heat and corrosion. It is fully insulated with glassfibre to avoid the overheating of external side. The door closes sterilization chamber hermetically by means of a silicon gum gasket resistant to heat. The sterilizer has five aluminium trays (four extractable trays and one as bottom surface) with holes for warm air circulation. The guides on which the trays run are separated from the sides to allow excellent warm air circulation (to obtain maximum heat distribution it is recommended not to use the bottom surface). There is a stainless steel electric immersion heater, characterized by high efficiency and little inertia. The automatic control of temperature is obtained with a capillary bulb automatic thermostat as auto-adjustment up to 200°C. If the temperature exceeds the limit set up during the test, a second safety bimetallic thermostat takes control. The working of the second thermostat causes the alarm red pilot lamp lighting up. To program the thermostat starting time, it is possible to use the timer up to 120 minutes or fixed on the manual position. The required temperature and time setting, is obtained by means the thermostat and timer knobs.

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
The device is provided with a green bright main power switch and with green, orange and red pilot lamps to show respectively the working, heating and alarm signals. The sterilizer is delivered with certified feeder cable, four extractable aluminium trays, two holder-tray pliers.

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

The sterilizer is shipped in an adequate package which requires the common handling cautions for fragile materials, as shown on the package itself. Use only the original package if the equipment must be returned to the manufacturer. To carry the sterilizer without package, keep it in a horizontal position and use the base as holding point.

Note : the cardboard package doesn't require particular waste disposal conditions.

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Remove the equipment from the package handling with care. Check the power supply system is correctly grounded. Place the sterilizer on a level surface so that air can freely circulate under the equipment bottom side. Keep the back side of the sterilizer at least 10cm. away from walls to ensure an adequate ventilation. Insert the trays in their housing. Insert the safety plug of the power supply cable, firstly into the equipment socket and then into the electric system, after being sure the voltage corresponds to the one selected on the device. Press the green main power switch position "I". Check if the switch light is on.

Note : in case you find a damp package, it is recommended to leave the unit in a heated room to exclude humidity conditions of the components.

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Note : *a typical smell of heated glassfibre may be present at the first power on and remain for a short period of time.*

Note : *you can find the environmental required conditions (pressure, humidity, temperatures, etc.) inside the technical data page, as the reference rules.*

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For a correct use and a good circulation of hot air, we recommend do not fill up the trays excessively, and we recommend do not put one instrument on another. The maximum weight that every tray can support is approx.: " Gimette 21" : 4,5 kg. – " Gimette 28" : 3 kg. – " Gimette 50" : 1,5 kg.

When the required temperature is reached, the sterilization times –according to the material- vary from min. 30 minutes up to max.120 minutes with temperature 150÷180°C max. according to the directions of "Farmacopea Ufficiale" and the medical devices manufacturers. The user must to value if the instruments' quality and characteristics, during the process, can exclude a deterioration of the products which are treated. It is important to specify the sterilization process is just only a part included in a more complex procedure which concerns the Decontamination, the Cleansing, the Sterilization and the Curing phases. The correct application of this procedure, the periodical activities of the chamber cleaning, can be exclude a crossed contamination. Do not keep the thermostat highest temperature for a long time, in order not to risk of spoiling the sterilized material. The lighting up of the red light (nr.3) "Alarm" indicates the temperature is going over 220°C and the second safety thermostat is working limiting the temperature. Let the service or the specialized technician verify the efficiency of the adjustable thermostat. About the sterilizers models without a thermometer it will be opportune a periodical check of the thermic performance to warrant the set up processes.

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Note :

periodically check with proper indicators, the heat homogeneous distribution inside the chamber.

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Always check the power is unplugged before maintenance is carried out. The sterilizer needs no particular maintenance, provided the adequate environment are satisfied, except for a regular cleaning of its external parts. It is recommended to clean the equipment with a soft cloth and water. Do not use acids or solvents. Regularly check if the power cable is in good conditions.

It is possible to replace the safety fuse pulling out the fuse block at the rear under the socket. Please always follow the characteristics stated in the technical data sheet. If the sterilizer doesn't work properly, contact the manufacturer or the authorized service. If any repairs are needed, please return the equipment directly to the manufacturer or to the local distributor using only the original package. However it is possible to contact specialized technicians, who must work in accordance to general safety and accident prevention rules of the country in which the equipment is installed.

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4.5 Messa fuori uso / *Putting out of service / Mise hors service*

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No particular cautions. In any case, about the waste process, it must be applied the required procedures by the rule in force and the EU Directives.

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Even though the sterilizer doesn't fall inside the table of the devices provided of the Annex RAEE Directive, in any case the unit has been conceived with a selection of components according to such Directive.

N.B.

From the technical point of view, the life of the device is related to the characteristics of the components used to manufacture the device itself, and it is not related to its execution from Fast Assembler side.

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This equipment must be only used from authorized personnel informed on the correct use.

These sterilizers must be only used for their designed purpose : to sterilizer surgical instruments or tools in the medical and paramedical field.

The sterilizers must not be used for all treatments concerning contact lenses (Rule 15 - All.9 - Comma 1).

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It is necessary to take all precautions to avoid unauthorized use. Do not modify, change or put out of order any electrical or thermal safety protection of the equipment. Check regularly the functionality of the light indicators which signal the feeding and the heating phases. When the sterilizer is operating, pay attention to the people near the equipment. The user must not attempt to carry out any repair unless authorized. Do not clean or repair the sterilizer while it is connected to an electrical system and be sure the sterilization chamber is at room temperature. Check regularly if the power cable is damaged. When authorized technicians decide to replace spare parts, they must use original replacements and moreover they must adhere to the manufacturer's information and technical characteristics. In case of any doubt about safety rules, please contact the manufacturer or local distributor before using the equipment. However use the sterilizer in accordance with the general safety and accident prevention rules of the country in which the equipment is installed.

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Avoid the accidental contact with liquids.
The sterilizer does not have to be used in an environment where the explosive gases may be present.

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Place the sterilizer on a level surface so that air can freely circulate under the equipment bottom side. Keep the sterilizer at least 10cm. away from walls to ensure an adequate ventilation. Keep the door closed when the sterilization process is running. Do not place any object on the upper part of the equipment. To extract the trays and fit them in, use the appropriate holder-tray pliers. Do not clean, move or repair the sterilizer while is connected to an electrical system and the sterilization chamber is still hot.

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When the sterilizer is operating do not put hands inside the sterilization chamber and do not touch the trays.

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The sterilizer must be used only for its designed purpose. Improper use of the equipment is to the risk and sole responsibility of the user.

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The warranty covers only equipment with production and natural defects as indicated in the certificate that comes with the device.

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