

PLANT

Upmarket podiatry chair equipped with 3 motors that control height, backrest and seat inclination. Reclining armrests and headrest cushion adjustable in multiple fixed positions. Legrests swivel in two fixed positions up to 90° and are independently adjustable in height by gas lift. White or grey PU upholstery of high quality that eases cleaning.

Sillón de podología de alta gama equipado con 3 motores que controlan la altura, la inclinación del respaldo y del asiento. Reposabrazos abatibles y cojín reposacabezas ajustable en múltiples posiciones prefijadas. Reposapiés ajustables en altura independientemente mediante pistón de gas y con apertura lateral en dos posiciones prefijadas hasta 90°. Tapizado blanco o gris en PU de alta calidad y fácil de limpiar.

Fauteuil de podologie haut de gamme équipé de 3 moteurs qui contrôlent la hauteur et l'inclinaison du dossier et du siège. Repose-bras rabattables et coussin repose-tête réglable en plusieurs positions fixes. Hauteur des jambières réglable de façon indépendante par vérin à gaz et écartables en deux positions jusqu'à 90°. Revêtement en PU blanc ou gris de haute qualité et facile à nettoyer.

Hochwertige Fußpflegeliege mit 3 Motoren, um die Höhe, die Neigung der Rückenlehne und die Neigung des Sitzes zu regulieren. Abklappbare Armlehnen und Kopfteil in festen Stellungen regulierbar. Fußstützen mittels Gasfeder unabhängig voneinander in der Höhe regulierbar und mechanisch um 90° drehbar. Hochwertiger weißer oder grauer PU-Polsterbezug, leicht abwaschbar.

Poltrona podologica di alta gamma dotata di 3 motori che controllano l'altezza, l'inclinazione dello schienale e del sedile. Braccioli pieghevoli e cuscino poggiatesta regolabile in più posizioni preimpostate. Poggiapiedi regolabili in altezza in modo indipendente tramite pistone a gas e con apertura laterale in due posizioni preimpostate fino a 90°. Rivestimento bianco o grigio in PU di alta qualità e facile da pulire.

Electric Podiatry Chair Ref. WKP005.3

Size: 168x81,5x54/130 cm
Weight: 79 Kg
Upholstery: PU
Motors: 3



- 1 Angle 0° ~ 80°
- 2 Extension ±18 cm
- 3 Angle 0° ~ 80°
- 4 Swivel legrest angle 0° ~ 90°
- 5 Angle 0° ~ 12°
- 6 Up and Down 54/128 cm
- 7 Extension ±10 cm

OPTIONS

- Bed cover 2232A (5 pieces)
- 8 Foot controller Ref. FC003.WKP005
- Wheels Ref. WH01.2
- Paper roll holder Ref. 50080.1
- 9 Debris tray Ref. 30030

COLOURS



PACKING SIZE cm

148x73x65

N.W. kg

77

G.W. kg

91

MAX WEIGHT

180 kg



PLANT

