

LUPA H.F

Cold light luminaire with magnifying lens, essential for high precision work. An ideal use for hospitals, medical centers, laboratories and aesthetic clinics.

Light Source

Low consumption 22W circular fluorescent tube.

Produces a light without dark spots and does not overheat.

Lifetime of 10.000 hours.

Functionality and performances

Adjustable height articulated arm that allows easy movement in different positions.

Head with biconvex circular lens of 5 diopters, adjustable and simple mechanism.

Switch on the head of the luminaire.

Security and hygiene

Head with protective cover to prevent dust or scratches that may damage the lens.

Mounting Options

Standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.



LUPA H.F

Technical Details	(Cod. 09332)
Light source	FL 22W
Lamp lifetime	10.000 hours
Colour temperature	6.500 °K
Glass lens	5 diopters
Lens diameter	125 mm
Arm length	90 cm
Weight	3 kg
Colour	White

Trolley stand
8,8 kg Cod. 09651





Wall support
Cod. B0101800



AH table clamp
Cod. B010190



Rail Plus bracket
Cod. 82026



Extension wall arm
Cod. B0101S0

