

# LUXIFLEX

Luxiflex is a versatile and stylish lamp for examination and minor surgery procedures, ideal in health centers, clinics and hospitals.

It is also suitable for veterinarian examinations and procedures.

Flexible arm allows to adjust the height to assure the proper working distance.

Low temperature lamp without infrared radiation. Switch located on the arm's base.

It can be adapted to: trolley stand, table clamp (included in standard version), rail support, wall support and extension wall arm.

Also available in LED version: Luxiflex LED and Luxiflex LED Plus.



## LUXIFLEX HALOGEN

Technical Details		(Cod. I0104AO)
Light source	Hal MasterLine ES 35W 8°	
Lamp lifetime	5.000 hours	
Lux a 50 centimeters	50.000 lux	
Colour temperature	4.000 °K	
Colour rendering	95%	
Total length	85 cm	
Colour	White	
Glass filter for colour temperature		



**Wall support**  
Cod. B0101800



**AH table clamp**  
Cod. B010190



**Rail Plus bracket**  
Cod. 82026



**Extension wall arm**  
Cod. B0101S0



**Trolley stand**  
4 kg Cod. 09655  
8,8 kg Cod. 09651