

# LUXIFLEX LED

Luxiflex LED is a cold light luminaire, versatile with a compact design. Luxiflex is designed for examination procedures that require great performance.

## Flexibility

Flexible arm allows to adjust the height to assure the proper working distance.

## Compact design

Luxiflex offers easy use, thanks to its compact LED head with a diameter of 8 cm and its switch located at the base of the lamp. These features make this model an ideal luminaire for the most demanding examinations procedures.

## Mounting options

The standard version includes AH table clamp.


As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.

## High performance LED

Luxiflex LED 15.000 lux at 50 centimeters.

## LUXIFLEX LED

Technical Details (Cod. I0104AOLED)	
Light source	LED 6W 15°
Lamp lifetime	50.000 hours
Lux a 50 centimeters	15.000 lux
Colour temperature	4.500 °K
Colour rendering	80%
Total length	85 cm
Colour	White





**Wall support**  
Cod. B0101800



**AH table clamp**  
Cod. B010190



**Rail Plus bracket**  
Cod. 82026



**Extension wall arm**  
Cod. B0101S0



**Trolley stand**  
4 kg Cod. 09655  
8,8 kg Cod. 09651