

BEDSIDE CABINETS



RUBENS 1

Practical design and a wide range of decors make the Rubens 1 bedside cabinet a perfect solution for storing personal belongings. It provides an excellent complement to the Elbur care beds. Convenient table top folding system, drawers and doors opened on both sides make that the use of the cabinet is very simple and comfortable.

Safe working load: 6 kg

Dimensions (height × depth × width) [cm]: 87×45×53 (119)

Weight: 34 kg



Rubens 1

RUBENS 7

Rubens 7 is a perfect complement to the Elbur care beds. This bedside cabinet has a drawer, an open shelf and a closed space with a shelf that allows patients to store their personal belonging within easy reach from the bed.

Dimensions (height × depth × width) [cm]: 87×45×43

Weight: 25 kg



Rubens 7

RUBENS 2

Rubens 2 bedside cabinet is distinguished by its great versatility thanks to additional height and angle adjustment system of the wooden table top. Invaluable functionality along with the possibility of a harmonious decor fit with the Elbur care beds, make this bedside cabinet model a very frequent choice as an equipment in the long-term care facilities and in home care as well.

Safe working load: 6 kg

Dimensions (height × depth × width) [cm]: 87×45×53 (121)

Weight: 41 kg

Height adjustment range [cm]: 71÷99

Angle adjustment range [°]: 0÷35



Rubens 2: the cabinet with height and angle adjustment system of the wooden side top

OVERBED TABLES



RUBENS 3

Rubens 3 overbed table works perfectly as a compliment to our beds. A wide range of height adjustment and multi-angle adjustable system of the wooden top simplify patients eating meals, reading or using mobile devices.

Safe working load:	11 kg
Dimensions (height × depth × width) [cm]:	72×90×45
Weight:	15 kg
Height adjustment range [cm]:	72÷107
Angle adjustment range [°]:	0÷52



Rubens 3



Rubens 9

RUBENS 9

An economical model of an overbed table is a proven solution in the long-term care. Angle adjustment of a wooden top and a wide range of height adjustment, improve user comfort and quality patient care. Rubens 9 stands out as a very stable and high quality product. It has 6-position height adjustment in the range of approx. 79-114 cm.

Safe working load:	20 kg
Dimensions (height × depth × width) [cm]::	79×78×42
Weight:	10 kg
Height adjustment range [cm]:	79÷114
Angle adjustment range [°]:	0÷45