

Work happily



V1.0

naggura'swop


naggura.com

Our products



Medical Examination

S301




Swop 301 is an all-round electrical treatment table! Simply multipurpose!

Wood, Metal or Pmma




Collections


Nordic
METAL



Yellow




Blue




Pink

Foll
METAL




Ash




Clay

Iron
METAL



Inox



Black

Om
WOOD



Oak




Wenge

Techno
PMMA



Fluor




Fuxia




Ice

S301



Dimensions



General Specifications

MOTORS

1 engine (height)

DIMENSIONS (LENGTH X WIDTH)

195 x 68 cm

HEIGHT (MIN. AND MAX.)

51/92 cm

WEIGHT

92 kg

THICKNESS

6 cm

DENSITY

30 kg/m³

LIFTING CAPACITY

175 kg

OPERATION

Hand controller (included)/
Foot controller (optional)

POWER SUPPLY

230 V - 50/60 Hz

MAX. INTENSITY

2 A

ELECTRIC MOTOR

6.000 N

CUSHION

Seamless Upholstery and
round edges.

100% Vinyl-compound

Permablok 3 advanced vinyl
protection (against germs,
abrasion), UV resistant, supe-
rior abrasion 300K Grey colour

Permaguard advanced vinyl
protection (against germs and
stains), UV resistant, abrasion
100K White colour