

TECHNICAL INFORMATION

FIX / FIX HL-LOW GYNAECOLOGICAL TABLE

Total dimensions	185 x 60 x 65/90 cm.
Back dimensions	90 x 60 cm.
Pelvis section dimensions	41 x 60 cm.
Back inclination	70°
Trendelenburg inclination	20°
Weight	62 kg.
Product manufactured in accordance with directive	93/42/CE
Class	1

NOTE: The measurements and characteristics of the products shown may be modified without notice.



 **medisa**
cares for you

Clinical Furniture
Manual gynaecological
examination tables



www.medical-iberica.com

Lanzahita, 6 - Pol. Ind. Albarreja - 28946 Fuenlabrada, Madrid (Spain) - T. 91 208 94 00 - F. 91 208 94 01

FIX examination table

A functional, economical solution at the service of gynaecologists

The FIX examination table has been designed for all medical surgeries that require a practical product in technical and economic terms.



STRUCTURE

- Structure made of steel tube with an epoxy coating.
- Base resting on 4 levelling legs.
- Fixed height.



Ref. MB-GINE-05
Leg and heel supports

COMFORT

- The Fix table has been submitted to stringent laboratory tests and trials in order to obtain the necessary degree of ergonomic adaptation and maximum comfort required to improve the quality of life and the health of users.

SAFETY

- Table manufactured in accordance with directive 93/42/EC.
- Class M1 upholstery (flame retardant and antibacterial).

FIX HL-LOW examination table

A functional, economical solution at the service of gynaecologists



STRUCTURE

- The base rests on 4 rotating wheels with brakes that have a diameter of 125 mm.
- Height adjustment system by means of a telescopic column.
- Minimum and maximum height: 83/123 cm.



HANDLING

- The back articulation and Trendelenburg position is activated by a gas piston located in the base which is easy and comfortable to use.
- Collapsible leg rest that are easy to remove and allow for a wide range of different configurations, such as horizontal positioning, in order to perform ultrasound examinations.

PADDING AND UPHOLSTERY

- Seat and back made of high density polyurethane foam, designed for optimum comfort and improved ergonomic adaptation Upholstered in self-extinguishing vinyl fabric (antibacterial, resistant to absorption and stains in accordance with the M1 standard). It is easy to clean and maintain.

ACCESSORIES

- Leg supports for patient with PUR padding. This material allows them to be adjusted in order to achieve a more comfortable position (included).
- Padded arm rests made of PUR, ref. MB-GINE-01.
- Roll holder ref. MB-GINE-02.
- Colposcopy support ref. MB-GINE-03.
- Adjustable proctology device ref. MB-GINE-04.