



NOTE: The measurements and characteristics of the products shown may be modified without notice.



Clinical Furniture

NOVA Model Modular Table



BUREAU VERITAS ECO-LABEL ISO 14001 IBV CE

www.medical-iberica.com

Nova Model Modular Table

Individual comfort and care thanks to its wide range of functions

Special treatments require special furniture.
Medisa provides functional, carefully designed products together with the highest quality.



DOUBLE MODULE TABLE

- Modular table with drawer unit and lectern
- Removable drawer unit with two compartments for improving the comfort of the patient while in bed.
- The drawer unit is designed to be hung on the wall.
- The table top can be used by the patient for different purposes such as eating or writing.
- The table is fitted with four wheels, two of which have a braking system.

Internal block	46 x 47 x 87,5 cm
External block	39 x 38,3 x 32,5 cm



DRAWER UNIT

- The drawer unit has separate wheels with a braking system.
- Removable single-piece drawer unit made of injection-moulded HD polyethylene that is easy to open and has a lock.

- Internal space with ventilation holes for improved hygiene.



TRAY (LECTERN)

- The height and angle of the tray are adjustable.
- It is laminated and has a non-scratch surface and raised edges for retaining liquids.
- Double compartment for bottles and a magazine and towel rack.

- Section for supporting diverse objects.
- Removable drawer with an internal mirror.

Tray	54 x 41 cm
Total measurements	78,5 x 59 cm

