

MIO-PERISTIM

Mio-Peristim is the electrostimulator designed to treat male and female incontinence problems.

It is able to stimulate the pelvic muscles, responsible for the control of the bladder, toning and strengthening them. It increases muscle control and decreases bladder irritability.

MEDICAL DEVICE CLASS IIa

MIO-PERISTIM is a medical device CE0068



KIT

- ✓ Mio-Peristim device
- ✓ TENS/NEMS/FES electrotherapy with graphic display
- ✓ 2 connection cables with splitters (8 electrodes connection)
- ✓ Pre-gelled adhesive electrodes
- ✓ Internal Ni-Mh rechargeable battery pack with new LL technology (long lasting)
- ✓ Battery charger
- ✓ Belt clip
- ✓ User manual and electrodes positions manual
- ✓ Carriage bag

TECHNICAL FEATURES

- ✓ Combined electrotherapy unit with 2 independent channels
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid hazardous thermal effect of polarization
- ✓ 9 INCONTINENCE (urinary and faecal) programs
- ✓ 14 TENS programs
- ✓ 13 free memories TENS - NEMS
- ✓ Digital easy to use keyboard
- ✓ Remaining therapy time and programs instructions shown on graphic display
- ✓ Internal rechargeable battery pack (removable)
- ✓ Backlight display
- ✓ Maximum intensity 200 mApp

DIMENSIONS

114 X 68 X 26 mm

MIO-PERISTIM

PROGRAMS

14 TENS PROGRAMS

Conventional tens (fast)

Endorphinic tens (delayed)

Tens at maximum values

Anti-inflammatory

Neck pain/cervicogenic headache

Backache/sciatic pain

Sprains/bruises

Vascularization

Muscle relaxant

Hand and wrist pain

Plantar stimulation

Epicondylitis

Epitroclea

Periarthritis

9 INCONTINENCE PROGRAMS

Stress urinary and faecal incontinence 1

Stress urinary incontinence 2

Stress urinary incontinence 3

Urge urinary and faecal incontinence 1

Urge urinary incontinence 2

Urge urinary incontinence 3

Mixed urinary and faecal incontinence 1

Mixed urinary incontinence 2

Mixed urinary incontinence 3

13 FREE MEMORIES

TENS mem 1

TENS mem 2

TENS mem 3

TENS mem 4

TENS mem 5

NEMS mem 1

NEMS mem 2

NEMS mem 3

NEMS mem 4

NEMS mem 5

NEMS alternated channels 1/2

NEMS alternated channels 1/2

Battery test

HOW DOES IT WORK

MIO-PERISTIM is a device particularly indicated for urinary and faecal incontinence treatments. It is provided with 9 incontinence programs using the vaginal and anal probe available on standard contents or as optional.

MIO-PERISTIM has also 14 TENS programs for pain treatment and 13 free memories to set TENS (antalgic therapy) and NEMS (neuromuscular stimulation) protocols.