

MS LED

MS LED is a minor surgery and advanced exploration luminaire. The lamp is manufactured with a LED light source and it is mounted with a highly precise articulated arm. Its arm allows the luminaire to be adjusted in any position.

The quality of the materials used in its manufacture makes it especially suitable for intensive use in surgery and examinations.

Exceptional quality of LED light

The powerful LED light source offers exceptional light quality and a lifespan of more than 50,000 hours.

The high performance of the colour Ra 97% depicts the exceptional quality of white light for an excellent distinction of fabrics and its reproducing all visible colours, naturally and precisely.

Ergonomic design

The circular shape of the head and its handle in the lower part facilitate its exact positioning, fulfilling the ergonomic requirements for any use.

Designed with clean and soft shapes which facilitate hygiene and cleanliness in medical environments.

Mounting options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.

MS LED

Technical Details	(Code M183500)
Light source	LED 8W
Lux a 50 centimeters	25.000 lux
Adjustable light intensity	No
Colour temperature	4.500 °K
Colour rendering	> 97%
Lamp lifetime	50.000 hours
Arm length	90 cm
Colour	White





Wall support
Code B0101800



AH table clamp
Code B010190



Rail Plus bracket
Code 82026



Extension wall arm
Code B0101S0



Trolley stand
8,8 kg Code 09651

