

ELBUR PB 325

(PB325L-4C-90200-I, article no.: 82200)

TECHNICAL SPECIFICATION		
Load [kg]	Max. Patient weight	140
	Safe Working Load	175
	Lifting pole load	80
Dimensions [cm]	External [W x L]	approx. 103 x 216
	Mattress support platform [W x L]	90 x 200
	Under bed clearance	> 15
Adjustment	Height [cm]	approx. 40,0 ÷ 80,0
	Angle – back rest section [°]	approx. 0 ÷ 70
	Angle – upper leg rest section [°]	approx. 0 ÷ 32
	Angle – lower leg rest section [°]	approx. 0 ÷ 20
Mattress	Dimensions [HxWxL]	10x90x200 / 12x90x200
	Volume density of foamed material	35 - 50 kg/m ³
Drive system	Drive lifting unit	Linak LA27
	Mattress support platform drive unit	Linak LA27
	Control box	Linak CB6
	Hand switch	HL74
	Input voltage	230V AC, 50Hz
	Operating time (duty cycle)	2 min on / 18 min off
	Ingress protection	IP X4
	Electrical appliance class	II
	Noise level	< 65 dB (A)
	Weight [kg]	Total
Mattress support platform (head section)		16,5
Mattress support platform (foot section)		16,3
End board with lifting drive unit (1)		15,7
Mattress support platform drive unit: back rest / leg rest		1,5
Hand switch / control box		0,3 / 1,6
Lifting pole (optional)		11,2
Conditions operating		Temperature [°C]
	Humidity [%]	+10 ÷ +40
Conditions storage	Temperature [°C]	30 ÷ 75
	Humidity [%]	+5 ÷ +45



According to the condition of **26th February 2021**

ELBUR PB 326

(PB326L-4C-90200-I, article no.: 82400)

TECHNICAL SPECIFICATION		
Load [kg]	Max. Patient weight	140
	Safe Working Load	175
	Lifting pole load	80
Dimensions [cm]	External [W x L]	approx. 104 x 216
	Mattress support platform [W x L]	90 x 200
	Under bed clearance	> 15
Adjustment	Height [cm]	approx. 40,0 ÷ 80,0
	Angle – back rest section [°]	approx. 0 ÷ 70
	Angle – upper leg rest section [°]	approx. 0 ÷ 32
	Angle – lower leg rest section [°]	approx. 0 ÷ 20
Mattress	Dimensions [HxWxL]	10x90x200 / 12x90x200
	Volume density of foamed material	35 - 50 kg/m ³
Drive system	Drive lifting unit	Linak LA27
	Mattress support platform drive unit	Linak LA27
	Control box	Linak CB6
	Hand switch	HL74
	Input voltage	230V AC, 50Hz
	Operating time (duty cycle)	2 min on / 18 min off
	Ingress protection	IP X4
	Electrical appliance class	II
	Noise level	< 65 dB (A)
	Weight [kg]	Total
Mattress support platform (head section)		17,0
Mattress support platform (foot section)		16,6
End board with lifting drive unit (1)		17,2
Mattress support platform drive unit: back rest / leg rest		1,5
Hand switch / control box		0,3 / 1,6
Lifting pole (optional)		11,2
Conditions operating		Temperature [°C]
	Humidity [%]	+10 ÷ +40
Conditions storage	Temperature [°C]	30 ÷ 75
	Humidity [%]	+5 ÷ +45



According to the condition of **26th February 2021**

ELBUR PB 331

(PB331L-4C-90200-I, article no.: 82600)

TECHNICAL SPECIFICATION		
<i>Load</i> [kg]	Max. Patient weight	140
	Safe Working Load	175
	Lifting pole load	80
<i>Dimensions</i> [cm]	External [W x L]	approx. 102 x 215
	Mattress support platform [W x L]	90 x 200
	Under bed clearance	> 15
<i>Adjustment</i>	Height [cm]	approx. 40,0 ÷ 80,0
	Angle – back rest section [°]	approx. 0 ÷ 70
	Angle – upper leg rest section [°]	approx. 0 ÷ 32
	Angle – lower leg rest section [°]	approx. 0 ÷ 20
<i>Mattress</i>	Dimensions [HxWxL]	10x90x200 / 12x90x200
	Volume density of foamed material	35 - 50 kg/m ³
<i>Drive system</i>	Drive lifting unit	Linak LA27
	Mattress support platform drive unit	Linak LA27
	Control box	Linak CB6
	Hand switch	HL74
	Input voltage	230V AC, 50Hz
	Operating time (duty cycle)	2 min on / 18 min off
	Ingress protection	IP X4
	Electrical appliance class	II
	Noise level	< 65 dB (A)
	<i>Weight</i> [kg]	Total
Mattress support platform (head section)		17,0
Mattress support platform (foot section)		16,6
End board with lifting drive unit (1)		24,9
Mattress support platform drive unit: back rest / leg rest		1,5
Hand switch / control box		0,3 / 1,6
Lifting pole (optional)		11,2
<i>Conditions operating</i>		Temperature [°C]
	Humidity [%]	+10 ÷ +40
<i>Conditions storage</i>	Temperature [°C]	30 ÷ 75
	Humidity [%]	+5 ÷ +45



According to the condition of **26th February 2021**

ELBUR PB 337

(PB337L-4C-90200-S, article no.: 83050)

TECHNICAL SPECIFICATION		
Load [kg]	Max. Patient weight	140
	Safe Working Load	175
	Lifting pole load	80
Dimensions [cm]	External [W x L]	approx. 102 x 215
	Mattress support platform [W x L]	90 x 200
	Under bed clearance	> 15
Adjustment	Height [cm]	approx. 23,0 ÷ 63,0
	Angle – back rest section [°]	approx. 0 ÷ 70
	Angle – upper leg rest section [°]	approx. 0 ÷ 32
	Angle – lower leg rest section [°]	approx. 0 ÷ 20
Mattress	Dimensions [HxWxL]	10x90x200 / 12x90x200
	Volume density of foamed material	35 - 50 kg/m ³
Drive system	Drive lifting unit	Linak LA27
	Mattress support platform drive unit	Linak LA27
	Control box	Linak CB6
	Hand switch	HL74
	Input voltage	230V AC, 50Hz
	Operating time (duty cycle)	2 min on / 18 min off
	Ingress protection	IP X4
	Electrical appliance class	II
	Noise level	< 65 dB (A)
	Weight [kg]	Total
Mattress support platform (head section)		20,3
Mattress support platform (foot section)		18,0
End board with lifting drive unit (1)		24,9
Mattress support platform drive unit: back rest / leg rest		1,5
Hand switch / control box		0,3 / 1,6
Lifting pole (optional)		11,2
Conditions operating		Temperature [°C]
	Humidity [%]	+10 ÷ +40
Conditions storage	Temperature [°C]	30 ÷ 75
	Humidity [%]	+5 ÷ +45



According to the condition of **26th February 2021**

ELBUR PB 331 & PB 337 with axional brakes

TECHNICAL SPECIFICATION		
Load [kg]	Max. Patient weight	140
	Safe Working Load	175
	Lifting pole load	80
Dimensions [cm]	External [W x L]	approx. 102 x 215
	Mattress support platform [W x L]	200
	Under bed clearance	> 15
Adjustment	Height [cm]	see Tab.2.
	Angle – back rest section [°]	approx. 0 ÷ 70
	Angle – upper leg rest section [°]	approx. 0 ÷ 32
	Angle – lower leg rest section [°]	approx. 0 ÷ 20
Mattress	Dimensions [HxWxL]	10x90x200 / 12x90x200
	Volume density of foamed material	35 - 50 kg/m ³
Drive system	Drive lifting unit	Linak LA27
	Mattress support platform drive unit	Linak LA27
	Control box	Linak CB6
	Hand switch	HL74
	Input voltage	230V AC, 50Hz
	Operating time (duty cycle)	2 min on / 18 min off
	Ingress protection	IP X4
	Electrical appliance class	II
	Noise level	< 65 dB (A)
	Weight [kg]	Total
Mattress support platform (head section)		20,3
Mattress support platform (foot section)		18,0
End board with lifting drive unit (1)		24,9
Mattress support platform drive unit: back rest / leg rest		1,5
Hand switch / control box		0,3 / 1,6
Lifting pole (optional)		11,2
Conditions operating	Temperature [°C]	5,3
	Humidity [%]	+10 ÷ +40
Conditions storage	Temperature [°C]	30 ÷ 75
	Humidity [%]	+5 ÷ +45



Tab.2. Height adjustment range

No.	Castors diameter	PB 331 with central brake	PB 337 with central brake
1.	φ10,0 cm	40,0 ÷ 80,0 cm	23,0 ÷ 63,0 cm
2.	φ 12,5 cm	44,0 ÷ 84,0 cm	27,0 ÷ 67,0 cm
3.	φ 15,0 cm	46,5 ÷ 86,5 cm	29,5 ÷ 69,5 cm

According to the condition of **26th February 2021**