



# NORSE RESCUE® EMERGENCY CHEST DRAINAGE KIT

In the most demanding medical emergencies, every second counts. The NORSE RESCUE® Emergency Chest Drainage Kit is meticulously designed to provide rapid, efficient, and safe chest drainage in cases of traumatic pneumothorax, hemothorax, and other life-threatening situations.

The kit contains the following items:

**Drainage Set:** A closed system featuring a CH28, 100 cm tube, and a 600 ml collapsible collection bag, ensuring a sterile environment for fluid collection.

**Safety Scalpel:** A size 22 blade with a safety cap, sterile, designed for precision and protection.

**Suture:** Polydox Monofilament, 150cm, violet, loop, with an HRX-48 needle (CTX), resorbable, and sterile for secure closure.

**Gloves:** Durable Nitrile gloves in size XL to ensure a safe and sterile working environment.

**Peen:** A 13 cm, straight, sharp, single-use, sterile clamp for controlling bleeding.

**Disinfection Wipes:** Large 140x190 mm wipes soaked in 82% Ethanol and 0.5% Chlorhexidine for thorough disinfection.

**Chest Seal:** A vented seal to prevent air from entering the pleural space while allowing air and fluid to escape.

