

# H 263

# hidemar

## ward screen with 3 rigid panels



long lasting



light weight



resistant



customizable



anti-static



antibacterial



humidity  
resistant



latex free

*Rigid panels made of phenolic resin in white colour, characterized by its high resistance to impacts, liquids and disinfectants; indestructible and easy to clean.*

*Independent folding panels.*

*ø50 mm double swivel casters.*

*Structure made of enameled steel epoxy coated (colour to choose)*

*Dimensions: 2100 x 400 x 1800 mm (length x width x height)*

*Panel surface: 650 x 1250 mm (width x height)*

*Weight: 22 kg*