

ELITE COMMODE WHEELCHAIRS - MADE IN ITALY



Same features as 43201 with adjustable height, without wheels



- 43201 ELITE 110 COMMUNE - 42 cm
- chromed steel frame
- blue, washable, flame retardant leatherette
- plastic pan
- 4 pivoting wheels Ø 8 cm (2 with brake)
- adjustable height footrest
- can be placed over WC

- 43203 ELITE 200 COMMUNE - 43 cm
- 43204 ELITE 200 COMMUNE - 53 cm
- chromed steel frame
- blue, washable, flame retardant leatherette
- plastic pan
- 4 wheels Ø 20 cm (2 with brake)
- adjustable length and raising footrest
- removable armrest
- reclining backrest

- 43199 ELITE 50 CHAIR
 - 43210 UNIVERSAL TRAY
 - 43212 WHEEL Ø 20 cm - front - spare
 - 43213 WHEEL Ø 20 cm - rear - spare
 - 43216 BELT
 - 43217 BELT with strap
- Sturdy and washable polypropylene belts. Quick mounting.

	TECHNICAL SPECIFICATIONS		
	ELITE 50	ELITE 110	ELITE 200
Size (cm):	52x60xh 70/78	52x80xh 100	54x88xh 105
Seat height	46/54 cm	50 cm	53 cm
Seat width	42 cm	42 cm	43 or 53 cm
Weight	20 kg	20 kg	21 kg
Load	90 kg	90 kg	100 kg

COMMUNE CHAIRS

- 27701 COMMUNE CHAIR - chromed steel
- Chromed steel chair frame, with plastic detachable commode seat.
- 27702 COMMUNE WHEELCHAIR - painted
- Powder coated steel frame, detachable armrest and footrest, solid rear castors with lock, plastic commode seat.
- 27703 COMMUNE WHEELCHAIR - stainless steel
- Stainless steel chair frame, detachable armrest and footrest, rear castors with lock, plastic commode seat.
- 27699 FOOTREST for 27702 - pair
- 27700 FOOTREST for 27703 - pair
- 27766 PLASTIC WC BOWL - spare - Ø 26xh 15 cm



GIMA code	Wheelchair icons						Commune icons		(Kg)	DIM	NW GW	
	55	44	12.5	12.5	51	88	90	44				42
27702/3	55	44	12.5	12.5	51	88	90	44	42	100	54x53x53	17 18
27701	56	44			47	94		43	38	100	50x19x75	8.5

WC - SHOWER COMMUNE



- 27758 BATH BENCH with backrest
- Aluminium frame and padded seat. Removable plastic bowl WC.
- 27766 PLASTIC WC BOWL - spare

- 27759 WC - SHOWER COMMUNE
- Suitable also as raised toilet seat
- Adjustable height 45-55 cm (back 78-88 cm)
- 27766 PLASTIC WC BOWL - spare

COMFORT COMMUNE CHAIR



The generously padded and removable seat offers this chair the best comfort and high level of hygiene. Pan is removable from both sides and from the top.

- 43190 COMFORT COMMUNE CHAIR - blue
 - 43191 COMFORT COMMUNE CHAIR - height adjustable - blue
- Made in France.



MAX LOAD 130 kg

TECHNICAL SPECIFICATIONS	
Total height: 86 or 86/97 cm	Height seat-floor: 52 or 52/63 cm
Width between armrests: 51.5 cm	Height seat-backrest: 38 cm
Total width: 60.5 cm	Seat plate depth: 39 cm
Height armrest-floor: 72 or 72/83 cm	Total depth: 58 cm