



1 pz.



**REF** 34328

**LOT** 0000xx

SET LARINGO "GIMA GREEN" F.O. LED 4 LAME MC-INT 1-2-3-4 - 2,5V  
 GIMA GREEN F.O.LED ADULT SET 4 BLADES MC-INTOSH 1-2-3-4 - 2.5 V  
 SET LED À F.O. GIMA GREEN ADULTE 4 lames: Mc Intosh 1-2-3-4 (2,5 V)  
 SET LARINGO "GIMA GREEN" F.O. LED ADULTOS 4 HOJAS MC-INT 1-2-3-4 (2,5 V)



GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy  
 Made in Pakistan - [www.gimaitaly.com](http://www.gimaitaly.com)



Size 19x23,5x4,8 cm  
509037 - Rev. 2 - 04/20

- REF** ○ **34328** Adult set 4 blades MClntosh 1-2-3-4 / Led bulb (2.5V)  
Set adult 4 blades MClntosh 1-2-3-4 / Lampadina Led (2.5V)
- REF** ○ **34329** Adult set 3 blades MClntosh 2-3-4 / Led bulb (2.5V)  
Set adult 3 blades MClntosh 2-3-4 / Lampadina Led (2.5V)
- Led Fiber Optic Laryngoscopes - Laringoscopi a Fibre Ottiche luce Led**

## Laringoscopi a Fibre Ottiche luce Led



- Laryngoscope à fibre optique - Led
- Laringoscópios em fibra ótica - Led
- Laryngoskopy światłowodowe - Led
- Λαρυγγοσκόπιο οπτικής ίνας - Led



- Luce fredda ad elevata luminosità - fibra Ø 4 mm
- Compatibile con tutti i sistemi F.O. "Green"
- Fibra facilmente rimovibile ed autoclavabile
- Luce Led

[www.gimaitaly.com](http://www.gimaitaly.com)



Gima S.p.A.  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)  
[www.gimaitaly.com](http://www.gimaitaly.com)  
Made in Pakistan

Vedere etichetta bar code  
See bar code label



**GIMA**  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)  
[www.gimaitaly.com](http://www.gimaitaly.com)

## Led Fiber Optic Laryngoscopes



- Laringoscopios de fibra óptica - Led
- Led LWL-Laryngoskope
- Laringoscop cu fibră optică - Led
- Led - مناظير الحجرة باللياف بصرية - Led



- Superior cold, bright light Ø 4 mm fibre
- Compatible with all "Green" F.O. systems
- Autoclavable, replaceable light guide
- Led light

[www.gimaitaly.com](http://www.gimaitaly.com)



**GIMA**  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)  
[www.gimaitaly.com](http://www.gimaitaly.com)