


Step-on personal floor scales KERN MPD · MPD-E

Professional step-on personal floor scales approved as a medical device and with EC type approval class III



KERN MPD 250K100M

Your KERN specialist dealer



Features

- Verification class III (verification is optional)
- Approved as a medical device according to 93/42/EEC
- The bestseller for doctor's surgeries, chemists and hospitals: easy to use with no waste of time – simply tap on it, stand on it, read the weight
- Robust design for daily use in a professional environment
- Ergonomically optimised design
- Easy and hygienic cleaning
- Secure and non-slip positioning with height-adjustable rubber feet
- Level indicator to level the balance precisely
- Large platform with non-slip and wear-resistant surface made of plastic, to provide a secure footing for patients
- Doctor scale MPD 200K-1EM: Identically to MPD 250K100M though without backlight and mains adapter, battery operation only, due to the reduced resolution faster display speed. Especially attractive in price!
- KERN MPD: Battery- or mains-powered, rechargeable battery operation optional
- Protective working cover included with delivery

Technical data

- Dimensions of weighing plate W×D×H 365×360×80 mm
- Dimensions of display device W×D×H 200×130×60 mm
- Net weight approx. 10 kg

MPD-E


- Large LCD display, without backlight, digit height 25 mm
- Battery operation, 6×1,5 V AA, operating time up to 20 h

MPD

- Large backlit LCD display, digit height 25 mm
- Battery operation possible, 6×1,5 V AA, not included, operating time up to 20 h
- Mains adapter external, standard

Accessories

- Protective working cover over the display device, scope of delivery: 5 items, KERN MBC-A06S05
- Cleaning cloths, alcohol-free cleaning cloths for wipe disinfection, quick-acting, on basis of modern quaternary ammonium compounds, and also against papova viruses. Particularly gentle on materials, and therefore very well suited for the disinfection of alcohol sensitive medical products. Fulfill the legal requirements for occupational safety in accordance with TRGS 525/540. Packaging contents 80 pcs., size 20×22 cm per cloth, KERN MYC-01
- KERN MPD: Rechargeable battery pack internal, operating time up to 40 h, without backlight, charging time approx. 12 h, KERN MBC-A08
- KERN MPD: Mains adapter external, 100-240 V, 50/60 Hz, standard EU, CH, UK, KERN YKA-43
- KERN MPD: Mains adapter external, 100-240 V, 50/60 Hz, standard EU, KERN YKA-44



* Within the EU, official verification is mandatory by law for scales that are intended for use as a medical device. Please add this to your order. We require the location of use and the post code for the verification.

Model	Weighing range [Max] kg	Readability [d] kg	Verification value [e] kg	Mandatory by law Verification	
				KERN	965-129
KERN					
MPD 250K100M	250	0,1	0,1	965-129	
MPD 200K-1EM	250	0,2	0,2	965-129	

Pictograms

Battery operation:
Ready for battery operation. The battery type is specified for each device

Rechargeable battery pack:
Rechargeable set

Universal mains adapter:
with universal input and optional input socket adapters for A) EU, CH

Verification possible:
The time required for verification is specified in the pictogram

Package shipment:
The time required for internal shipping preparations is shown in days in the pictogram