



Width
70cm



CLINNOVA Clinical 3 Section



These versatile 3 Section Couches provide a very stable support surface for prone, supine and seated treatment procedures and examinations. Both the elevating sections can be set to a wide range of angles for specific patient positioning, and are available in either electric or gas assisted.

Models

See Page 2 for all available models



Model NV3570-PRM
Colour Burgundy

Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Choice of hydraulic or electric height elevation
- Foot Section adjustment from 0° to +85°
- Head Section adjustment from -20° to +80°
- Electric models raise from low to high in only 15 seconds
- Excellent height range from 46cm to 95cm
- Linked retractable wheel design with adjustable foot
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

Optional Accessories

- 6004** Breathing hole & plug
- 9001** Perimeter foot switch operation
- 6200** Single foot switch for NV3560
- 6209** Triple foot switch for NV3570
- 9003** Side support rails (Pair)
- 6003P** Side support rail cushions (Pair)
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility



Low patient transfer height
Shown on NV3560-CLS
Colour Light Grey

Models Available



NV3550-CLS

Hydraulic Elevation

Gas Assisted Head & Foot Sections

Classic Base (CLS)



NV3560-CLS

Electric Elevation

Gas Assisted Head & Foot Sections

Classic Base (CLS)



NV3570-CLS

Electric Elevation

Electric Head & Foot Sections

Classic Base (CLS)



NV3550-PRM

Hydraulic Elevation

Gas Assisted Head & Foot Sections

Premium Base (PRM)



NV3560-PRM

Electric Elevation

Gas Assisted Head & Foot Sections

Premium Base (PRM)



NV3570-PRM

Electric Elevation

Electric Head & Foot Sections

Premium Base (PRM)