

Width
70cm



46cm to 95cm



CLINNOVA Mobile 2 Section



Very strong and rigid 2 section mobile treatment couches supplied with side support rails and a paper roll holder as standard. The Large Mobile Wheel Design incorporates central locking & steering facility and the couches have excellent height ranges with a low patient transfer height.

Models

See Page 2 for all available models



Model NV2660-CLS shown with optional side support cushions. Colour dark grey

Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Gas-assisted backrest adjustable from -20° to +80°
- Electric backrest feature on NV2470 & NV2670 (CLS & PRM)
- Choice of hydraulic or electric height elevation
- Electric models raise from low to high in only 15 seconds
- Excellent height range from 46cm to 95cm
- Choice of Individually Braking castors (IBC) or Large Mobile Wheel Design with central locking & steering facility (LMWD)
- Standard with side support rails and paper roll holder
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

Optional Accessories

- 6004** Breathing hole & plug
- 6200** Single Foot Switch for NV2460 & NV2660
- 6207** Double Foot Switch for NV2470 & NV2670
- 6002P** Cushions for side support rails (Pair)
- 6007** Removable head/knee support
- 6054** Telescopic IV pole

Model NV2450-PRM
Colour Sky Blue



Models Available



NV2450-CLS

Hydraulic Elevation
Gas Assisted Backrest
Independently Braking Castors
Classic Base (CLS)



NV2450-PRM

Hydraulic Elevation
Gas Assisted Backrest
Independently Braking Castors
Premium Base (PRM)



NV2460-CLS

Electric Elevation
Gas Assisted Backrest
Independently Braking Castors
Classic Base (CLS)



NV2460-PRM

Electric Elevation
Gas Assisted Backrest
Independently Braking Castors
Premium Base (PRM)



NV2470-CLS

Electric Elevation
Electric Backrest
Independently Braking Castors
Classic Base (CLS)



NV2470-PRM

Electric Elevation
Electric Backrest
Independently Braking Castors
Premium Base (PRM)



NV2650-CLS

Hydraulic Elevation
Gas Assisted Backrest
Central Locking Castors
Classic Base (CLS)



NV2650-PRM

Hydraulic Elevation
Gas Assisted Backrest
Central Locking Castors
Premium Base (PRM)



NV2660-CLS

Electric Elevation
Gas Assisted Backrest
Central Locking Castors
Classic Base (CLS)



NV2660-PRM

Electric Elevation
Gas Assisted Backrest
Central Locking Castors
Premium Base (PRM)



NV2670-CLS

Electric Elevation
Electric Backrest
Central Locking Castors
Classic Base (CLS)



NV2670-PRM

Electric Elevation
Electric Backrest
Central Locking Castors
Premium Base (PRM)