


Width
70cm



130cm

52cm

15cm



46cm to 95cm



CLINNOVA Therapy 2 Section – Plus



The plus head design is ideal for neck mobilisation as well as manual traction techniques. These couches feature a short head section with height adjustable armrests for prone treatments, promoting relaxed shoulders and a more comfortable patient. The couch has an excellent height range of 46 to 95cm and is fitted with a breathing hole & plug as standard.

Models

See Page 2 for all available models

Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Choice of hydraulic or electric height elevation
- Gas-assisted head section adjustment from -25° to +45°
- Head section features two individually adjustable armrests
- Electric models raise from low to high in only 15 seconds
- Excellent height range from 46cm to 95cm
- Linked retractable wheel design with adjustable foot
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

Optional Accessories

- 6200** Single Foot Switch
- 9002** Side Support Rails
- 6002P** Side Support Rail Cushions
- 9005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 9001** Electric Perimeter Foot Switch (EPFS)



Model NV2559-CLS
Colour Dark Grey



Low patient transfer height



Offset gas strut for unobstructed patient view



Models Available



NV2559-CLS

Hydraulic Elevation

Classic Base (CLS)



NV2569-CLS

Electric Elevation

Hand or Foot Switch

Classic Base (CLS)



NV2569-CLS-EPFS

Electric Elevation

Perimeter Foot Switch

Classic Base (CLS)



NV2559-PRM

Hydraulic Elevation

Premium Base (PRM)



NV2569-PRM

Electric Elevation

Hand or Foot Switch

Premium Base (PRM)



NV2569-PRM-EPFS

Electric Elevation

Perimeter Foot Switch

Premium Base (PRM)