



## CLINNOVA Therapy 2 Section – Standard



A very strong and rigid 2 section couch with a wide range of optional accessories available and an extensive choice of colours. This couch is suitable for supporting very heavy patients for examinations, minor procedures and physiotherapy treatments. This 2 section couch is available as a hydraulic or electric model, both models feature a gas assisted backrest that can be adjusted from -20° to +80° and has an excellent height range from 46cm to 95cm.

### Models

See Page 2 for all available models



Model NV2557-CLS  
Colour dark grey

### Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Choice of hydraulic or electric height elevation
- Gas-assisted backrest adjustment from -20° to +80°
- Electric models raise from low to high in only 15 seconds
- Excellent height range from 49cm to 99cm
- Linked retractable wheel design with adjustable foot
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

### Optional Accessories

- 6200** Single Foot Switch
- 9002** Side Support Rails
- 6002P** Side Support Rail Cushions
- 9005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 9001** Electric Perimeter Foot Switch (EPFS)

Negative head section to -20°  
Shown with optional paper roll holder



Offset gas strut for  
unobstructed patient view

Low patient transfer height  
Shown with optional paper roll holder



## Models Available

---



### **NV2557-CLS**

Hydraulic Elevation

Classic Base (CLS)



### **NV2567-CLS**

Electric Elevation

Hand or Foot Switch

Classic Base (CLS)



### **NV2567-CLS-EPFS**

Electric Elevation

Perimeter Foot Switch

Classic Base (CLS)



### **NV2557-PRM**

Hydraulic Elevation

Premium Base (PRM)



### **NV2567-PRM**

Electric Elevation

Hand or Foot Switch

Premium Base (PRM)



### **NV2567-PRM-EPFS**

Electric Elevation

Perimeter Foot Switch

Premium Base (PRM)