


MC7212
45cm to 59cm

MC7213
54cm to 74cm



Contoured Medical Chair



A versatile chair with a contoured seat pad and raised front to promote positive posture whilst seated. The 3 lever mechanism provides independent positioning of the backrest angle, seat pad tilt and chair height. The contoured backrest is ergonomically designed to provide lumbar support for the operator and has the option of foot height control.

Models

- MC7212** Contoured medical chair - standard model
- MC7213** Contoured medical chair - high model
Foot support ring fitted as standard



Model MC7212 Standard contoured medical chair



Model MC7213 High model contoured medical chair with **optional foot height control**