

## Economy Tilt Tables




These entry-level tilt tables have been designed for a wide range of rehabilitation, assessment and examination procedures. The table features 125mm individually braking castors, dual adjustable foot boards and a set of harness straps. For increased specification & features, please see the Tilt Table Pro.

## Models

- ST7640** Economy tilt table
- ST7640DL** Economy tilt table (Divided leg)



Individually adjustable dual foot boards



Model ST7640.  
Colour sky blue

## Standard Features

- Safe working load and lifting capacity of 225Kg (495lbs)
- 125mm individually braking castors
- Electric height and tilting operation with hand switch control
- Adjustable angle dual foot boards +15° to -30° (4 Positions)
- Height range of 47cm to 101cm
- Tilting angle from horizontal to +90°
- Lowers to wheelchair height for ease of patient transfer
- Tilt angle inclinometer fitted as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Supplied with patient handgrips and set of harnesses
- 70cm wide upholstery as standard

## Optional Accessories

- 6057** Battery Back-Up
- 6058** Lateral Supports (Pair)
- 6059** Adjustable Pommel Support\*
- 6007** Head Support Cushion
- 6062** Standard Harness Set
- 6063** Dual Leg Harness Set\*
- 6064** Pair Dual Knee Harnesses\*
- 6065** Standard Harness (Wide)
- 6066** Standard Harness (Narrow)
- 7648** Adjustable Work Table

\* (will fit ST7640DL only).



Angle inclinometer



Individually braked castors.  
Models ST7640 & ST7640DL



Model ST7640 in horizontal position.  
Colour sky blue