



Medicare 2 Section Couches



A very strong and rigid couch available with a wide range of accessories. Suitable for supporting very heavy patients for examination and minor treatment procedures.

Models

- SM2550** Hydraulic with manual backrest
- SM2560** Electric with manual backrest
- SM2570** Electric with electric backrest



Model SM2560 with optional foot switch.
Colour light grey

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable backrest from -25° to 75° on Model SM2550 and on Model SM2560
- Adjustable backrest from 0° to 85° on Model SM2570
- The power assisted backrest section incorporates anti-trap design with Models SM2550/2560
- Choice of manual or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design and adjustable foot for stability
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- Standard width 65cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6053** Protective base frame cover
- 6004** Breathing hole and plug for prone treatments
- 6001** Perimeter foot switch operation
- 6200** Single foot switch operation for Model SM2560
- 6207** Double foot switch operation for Model SM2570
- 6005** Paper roll holder
- 6002** Side support rails (Pair)
- 6002P** Side support rail cushions (Pair)
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole

See page 50-51 for additional information on our range of optional accessories



Low patient transfer height



Negative head section to -25°