



Medicare Bariatric 3 Section Couch



The Bariatric 3 Section Couch features electric height, backrest and foot section operation for effortless patient positioning for a wide range of examination procedures. There is a choice of wheel system available and the couch provides a safe working load and lifting capacity of 325Kg (715lbs).

Models

- SM3573** Electric with electric backrest and foot section. Retractable wheel system
- SM3673** Electric with electric backrest and foot section. Central locking wheel system



Model SM3573 Colour light grey

Standard Features

- Safe working load and lifting capacity of 325Kg (715lbs)
- Choice of retractable wheel or large wheel design, with central locking & steering facility
- Adjustable angle electric footrest operation from -25° to +45°
- Adjustable angle electric backrest from horizontal to +80°
- Electric height elevation from 45cm to 98cm
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Standard width upholstery of 80cm (maximum)
- Hand switch fitted as standard
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6076** Memory foam upholstery (60mm) for additional patient comfort
- 6004** Breathing hole and plug for prone treatments
- 6007** Head support cushion
- 6074** Side support rails (Pair)
- 6074P** Side support rail cushions (Pair)
- 6005** Paper roll holder

See page 48-49 for additional information on our range of optional accessories



Model SM3573 with optional breathing hole and plug. Colour light grey



Low patient transfer height