



Medicare Deluxe Gynaecology Couch



Designed for colposcopy and general gynaecology procedures, these Deluxe models feature electric positioning for height, backrest angle and tilt position. They are supplied with a special set of accessories which include a matching operator's chair, base frame cover, a choice of leg extension, leg supports, matching head support cushion, side support loops, debris tray and paper roll holder.

Models

- SM8583D** Deluxe Model with removable extension
- SM8593D** Deluxe Model with fold-down extension



Model SM8583D Deluxe model with removable extension. Colour cream

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable angle electric backrest from horizontal to +80°
- Electric control of height, backrest and tilt angle
- Excellent height range from 52cm to 106cm
- Crescent shaped seat design
- Fully adjustable patient leg supports
- Fitted with debris tray, side support loops and paper roll holder
- Supplied with matching operators chair and head cushion
- Supplied with the choice of removable or fold-down leg extension
- Compact base frame with retractable wheel system
- Base frame cover
- Adjustable foot for stability on uneven floors
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Hand switch fitted as standard
- Standard width 65cm upholstery
- Choice of 15 standard upholstery colours

Optional Accessories

- 6209** Triple foot switch
- 6057** Battery back-up facility



Fold-down leg extension



Model SM8593D Colour sky blue