



## Medicare Echocardiography Couch



Designed specifically for Echocardiography, this couch has electric operation of height, tilt and backrest angle. The backrest has a cut-away design for improved application of the transducer during scanning procedures and a specially designed stool for the clinician is also provided. The couch is also supplied with an integral footrest for the clinician, as well as a patient hand-grip and head cushion for the patients comfort and security.

### Models

- SM3590** Electric height, tilt and backrest. Retractable wheel system
- SM3690** Electric height, tilt and backrest. Central locking wheel system with steering facility



Model SM3590 Colour sky blue

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric adjustment of height, tilt and backrest angle with hand switch control
- Excellent height range from 50cm to 103cm
- Electric backrest angle adjustable from horizontal to +80°
- Electric trendelenburg tilt to -15°
- Manual adjustment of foot section from -25° to +70
- The foot section incorporates anti-trap design gas strut
- Backrest design with cut-away for improved transducer positioning
- Integral footrest for clinician when traditional right hand scanning, together with specially designed stool
- Hand-grip and side support with cushion for patient positioning and security
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Standard width 65cm upholstery
- Paper roll holder and head support cushion supplied as standard
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Choice of upholstery colours, for full range of colour options see page 51

### Optional Accessories

- 6209** Triple foot switch operation
- 6057** Battery back-up facility
- 6123** Fully adjustable arm support for the clinician

See page 48-49 for additional information on our range of optional accessories.



Special design stool provided for clinician, included as standard



Model SM3590 Colour sky blue



Model 6123 optional adjustable arm support for clinician