



Medicare Gynaecology Couches



Designed for general gynaecology procedures, these models are available with either hydraulic or electric height operation and feature a wide choice of optional accessories to meet both financial and clinical requirements. All models feature a crescent shaped seat design and are supplied with fully adjustable leg supports, paper roll holder, side support loops and debris tray for gynaecology procedures.

Models

- SM8553** Hydraulic with manual backrest
- SM8563** Electric with manual backrest
- SM8573** Electric height and backrest
- SM8583** Electric height, tilt and backrest



Model 6043 optional removable leg extension accessory



Model 6051 optional fold-down leg extension accessory



Model SM8583 Colour dark blue



Model SM8563 Colour light grey



Model 6100 see choice of matching operators chairs on page 50

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable angle electric backrest from horizontal to +80°
- Crescent shaped seat design
- Fully adjustable patient leg supports
- Fitted with debris tray, side support loops
- Fitted with paper roll holder
- Choice of manual or electric height elevation
- Excellent height range from 52cm to 106cm
- Compact base with retractable wheel system
- Adjustable foot for stability on uneven floors
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Hand switch fitted as standard on electric models
- Electric trendelenburg tilt to -15° on model SM8583
- Standard width 65cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6052** Protective base frame cover
- 6200** Single foot switch operation for Model SM8563
- 6207** Double foot switch operation for Model SM8573
- 6209** Triple foot switch operation for Model SM8583
- 6051** Fold-down leg extension accessory
- 6043** Removable leg extension accessory
- 6007** Matching head cushion (Included on Model SM8583)
- 6100** Matching operators chair
- 6057** Battery back-up facility

See page 48-49 for additional information on our range of optional accessories.