



## Medicare Multi-Couch Dual Footrest



Incorporating both an adjustable angle backrest with patient armrests and dual adjustable footrests, these models are suitable for a wide range of seated examination procedures such as podiatry, leg ulcer clinics, orthopaedic casting (plaster room), as well as being ideal for GP surgeries and health centres.

### Models

- SM0575** Electric height and backrest
- SM0585** Electric height, tilt and backrest



Model 6004 optional breathing hole and plug for prone treatments



Matching headrest as standard on Model SM0585. Option on Model SM0575



Model SM0585 with optional base cover. Colour sky blue



Low patient transfer height



Model 6150 optional matching premium operators chair

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric adjustable backrest from horizontal to +80°
- Adjustable angle footrests from horizontal to -85° and incorporate anti-trap design gas struts
- The footrests extend for tall patients
- The fold-down armrests have angular adjustment
- Electric height elevation via hand control
- Excellent height range from 52cm to 106cm
- Triple electric Model SM0585 includes paper roll holder and matching headrest as standard
- Electric tilt facility of -15° on Model SM0585
- Compact frame design with retractable wheel system
- Adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Standard width 65cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

### Optional Accessories

- 6052** Protective base frame cover
- 6007** Head support cushion (Included on Model SM0585)
- 6100** Matching operators chair
- 6150** Matching premium operators chair
- 6005** Paper roll holder (Included on Model SM0585)
- 6004** Breathing hole and plug for prone treatments
- 6209** Triple foot switch operation for Model SM0585
- 6207** Double foot switch operation for Model SM0575
- 6057** Battery back-up facility

See page 48-49 for additional information on our range of optional accessories.