



Medicare Multi-Couch Single Footrest



Incorporating both an adjustable angle backrest and a single drop end foot section, these models are suitable for a wide range of seated and supine patient examination procedures. The slightly narrow foot section enables the addition of leg supports or stirrups for gynaecological procedures.

Models

- SM9556** Hydraulic with manual backrest
- SM9566** Electric with manual backrest
- SM9576** Electric with electric backrest



Model SM9566 showing foot section at -85°. Colour light grey

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable angle backrest from horizontal to +80°
- Adjustable angle foot section from horizontal to -85°
- The power assisted backrest and foot sections incorporate anti-trap design gas struts
- Choice of manual hydraulic or electric height elevation
- Excellent height range from 50cm to 105cm
- Compact frame design with retractable wheel system
- Adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6004** Breathing hole and plug for prone treatments
- 6005** Paper roll holder
- 6006** Upper limb support
- 6007** Removable head/knee cushion
- 6052** Protective base frame cover
- 6054** Telescopic IV pole
- 6044** Lithotomy stirrups
- 6045** Leg supports for gynaecology
- 6057** Battery back-up facility
- 6200** Single foot switch operation for Model SM9566
- 6207** Double foot switch operation for Model SM9576

See page 48-49 for additional information on our range of optional accessories.



Low patient transfer height



Model 6045 optional leg supports for gynaecology



Model 6044 optional lithotomy stirrups