



Medicare Orthopaedic Couches



Designed specifically for plaster casting technicians, this couch features individually adjustable leg sections and is available with a range of optional accessories to assist and support the technician in carrying out various casting procedures. The Deluxe version is supplied as a package with base shroud, set of accessories and matching stool.

Models

- SM5555** Hydraulic with manual backrest
- SM5565** Electric with manual backrest
- SM5575** Electric with electric backrest
- SM5585** Electric with electric height, tilt and backrest
- SM5585D** Deluxe model with stool, base cover and Orthopaedic accessories set



Model 5007 optional orthopaedic full leg seat module



Model 5006 optional orthopaedic arm support attachment



Model SM5585D
Colour black

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable angle backrest from horizontal to +80°
- Adjustable angle footrests from horizontal to -85°
- Trendelenburg tilt to -15° on Models SM5585/ SM5585D
- The power assisted foot sections incorporate anti-trap design gas struts
- The footrests extend for tall patients
- Choice of manual hydraulic or electric height elevation
- Excellent height range from 52cm to 106cm
- Compact frame design with retractable wheel system
- Adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6052** Protective base frame cover
- 5006** Orthopaedic arm support
- 5007** Orthopaedic full leg seat module
- 5008** Knee support roll set
- 5009** Orthopaedic leg support attachment
- 6007** Removable head/knee cushion
- 6103** Matching operators stool
- 6107** Deluxe saddle chair
- 6005** Paper roll holder
- 6004** Breathing hole and plug for prone treatments
- 6200** Single foot switch operation for Model SM5565
- 6207** Double foot switch operation for Model SM5575
- 6209** Triple foot switch operation for Model SM5585 / SM5585D
- 6057** Battery back-up facility

See page 48-49 for additional information on our range of optional accessories.



Model 5008 optional knee support roll set



Model 5009 optional orthopaedic leg support attachment



Model SM5565 showing extended patient footrest