



Height range 52cm to 106cm



Height range 52cm to 106cm

## Medicare Phlebotomy Couches



A specialist couch designed to provide a comfortable and easy to mount seat position for blood sampling procedures. All models are fitted with a pair of multi-adjustable arm rests and side support loops for the stability and safety of the patient when mounting and dismounting the couch. The seat section is inclined for additional patient security.

### Models

<b>SM9556P</b>	Hydraulic with manual backrest
<b>SM9566P</b>	Electric with manual backrest
<b>SM9576P</b>	Electric with electric backrest
<b>SM9586P</b>	Electric height, tilt and backrest



Model SM9566P with optional base cover. Colour dark blue

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable angle backrest from horizontal to +80°
- Electric backrest and trendelenburg tilt on Model SM9586P. Electric backrest on Model SM9576P
- Patient hand-grip supports to aid mount/dismount from the couch as well as providing patient security
- Fully adjustable phlebotomy arms for use on either side of couch
- Fitted with paper roll holder
- Choice of hydraulic or electric height elevation
- Excellent height range from 52cm to 106cm
- Compact frame design with retractable wheel system
- Adjustable foot for stability
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Inclined seat section for patient security
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

### Optional Accessories

- 6052** Protective base frame cover
- 6200** Single foot switch operation for Model SM9566P
- 6207** Double foot switch operation for Model SM9576P
- 6209** Triple foot switch operation for Model SM9586P
- 6057** Battery back-up facility
- 6007** Removable head cushion (Included with model SM9586P)
- 6101** Matching high operators chair
- 6107** Matching deluxe saddle chair

See page 48-49 for additional information on our range of optional accessories.



Model 6107 optional matching deluxe saddle chair



Model 6101 optional matching high operators chair



Model SM9566P Colour dark blue