



Medicare Ultrasound Couches



Providing smooth electric positioning of height, tilt and backrest angle, these couches feature a compact base frame and are ideal for both supine and seated procedures. The deluxe model includes the base frame cover and matching operators chair. The new SM9696 models additionally feature an electric footrest for patient positioning. The slightly narrow foot section enables the additional fitting of leg supports or stirrups for gynaecological procedures.

Models

- SM9686** Electric height, tilt and backrest - manual foot section
- SM9686D** Deluxe Model with base cover and matching operators chair
- SM9696** Electric height, tilt, backrest and foot section
- SM9696D** Deluxe Model with base cover and matching operators chair

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Easy adjustment of height, tilt, backrest and foot section
- Electric adjustable backrest from horizontal to +80°
- Electric adjustable backrest from horizontal to +80°
- Adjustable angle footrest from horizontal to -85°
- Electric trendelenburg angle to -15°
- Choice of power assisted or electric operated foot section
- Electric height elevation via hand control
- Excellent height range from 52cm to 106cm
- Supplied with paper roll holder and matching head cushion
- Compact frame design with retractable wheel system
- Adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Standard width 70cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6052** Protective base frame cover
- MC6158** Matching sonographers chair
- MC6107** Matching deluxe saddle chair
- 6100** Matching operators chair
- 6004** Breathing hole and plug for prone treatments
- 6209** Triple foot switch for Model SM9686 & SM9686D
- 6057** Battery back-up facility
- 6045** Leg supports for gynaecology procedures
- 6044** Lithotomy stirrups
- 6054** Telescopic IV pole
- 6041** Side support loops

See page 48-49 for additional information on our range of optional accessories.



Model 6044 optional lithotomy stirrups



Model SM9686 with optional leg supports. Colour light grey



Model MC6158 Sonographers chair standard model. Colour oriflamme



Model SM9686 in trendelenburg position