



## Multi-Purpose Tilt Table



These versatile tilt tables feature a 2 section top, incorporating an adjustable backrest to meet a wide range of applications. Available with the option of "divided leg design" both models feature a large mobile wheel design with steering facility and have the standard benefits of dual adjustable foot boards, patient handgrips, adjustable worktable and a full set of harnesses. There is also a wide choice of optional accessories available to choose from.

### Models

- ST7647** Multi-Purpose Tilt Table with Gas Assisted Backrest
- ST7647DL** Multi-Purpose Tilt Table with Gas Assisted Backrest (Divided Leg)
- ST7649** Multi-Purpose Tilt Table with Electric Backrest
- ST7649DL** Multi-Purpose Tilt Table with Electric Backrest (Divided Leg)



Model ST7647. Colour sky blue

### Standard Features

- Safe working load and lifting capacity of 225Kg (495lbs)
- Large wheel design, with central locking & steering facility
- Adjustable angle backrest facility with anti trap design 0° to +75°
- Electric height and tilting operation with hand switch control
- Adjustable angle dual foot boards +15° to -30° (4 Positions)
- Height range of 47cm to 101cm
- Tilt table angle ranges from horizontal to 90°
- Lowers to wheelchair height for ease of patient transfer
- Tilt angle inclinometer fitted as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Supplied with patient worktable, handgrips and set of harnesses

### Optional Accessories

- 6057** Battery Back-up
- 6058** Lateral Supports (Pair)
- 6059** Adjustable Pommel Support\*
- 6007** Head Support Cushion
- 6062** Standard Harness Set
- 6063** Dual Leg Harness Set\*
- 6064** Pair Dual Knee Harnesses\*
- 6065** Standard Harness (Wide)
- 6066** Standard Harness (Narrow)

\* (will fit Divided Leg version only)



Angle inclinometer



ST7647DL Divided leg model



Model ST7647 in horizontal position with backrest adjusted to +75°