

**MC6112**  
45cm to 59cm



**MC6111**  
45cm to 59cm



**MC6113**  
45cm to 59cm



**MC6115**  
54cm to 74cm

**MC6114**  
54cm to 74cm

**MC6116**  
54cm to 74cm



## Laboratory Chairs & Stools



Designed for use within laboratories or clean room environments. Backrests and seat pads are all manufactured using hard wearing chemical resistant polyurethane. All chairs and stools are height adjustable and can be fitted with optional foot height control and a wide choice of castors.

### Models

- MC6111** Lab chair - standard model
- MC6112** Lab stool - standard model
- MC6113** Lab saddle stool - standard model
- MC6114** Lab chair with foot support ring
- MC6115** Lab stool with foot support ring
- MC6116** Lab saddle stool with foot support ring

### Features

- Designed for laboratories and cold rooms
- Hard wearing chemical resistant polyurethane (PU)
- Unique design improves ergonomic posture
- Wide choice of models available
- Restraint castors fitted on all models as standard

### Optional Accessories

- 6178** Anti-static castors
- 6455** Load locking castors
- 6172** Manual locking castors
- 6173** Glides for added stability on high models
- 6175** Height adjustable armrests
- 6174** Foot height control operation
- 6176** Foot support ring for standard model



Model MC6111 Laboratory chair, standard model



Model MC6114 Laboratory chair, high model with foot support ring



Model MC6113 Laboratory saddle stool



Model MC6115 Laboratory stool, high model with foot support ring