



Standard Tilt Tables



These excellent Tilt Tables have been designed to meet the requirements of a wide range of applications for rehabilitation, assessment and examination procedures. Available with the option of "divided leg design" which provides individual harness security for the lower limbs, and options of pommel and leg stabilisers. These models feature central locking castors with steering facility and include the adjustable work table, which is an option on economy models. All models are supplied with individually adjustable foot boards, patient handgrips and harnesses.

Models

- ST7641** Standard tilt table
- ST7641DL** Standard tilt table (Divided leg)



Individually adjustable dual foot boards
(Shown on Economy model)



Model ST7641.
Colour sky blue

Standard Features

- Safe working load and lifting capacity of 225Kg (495lbs)
- Large mobile wheel design, with central locking and steering facility
- Electric height and tilting operation with hand switch control
- Adjustable angle dual foot boards +15° to -30° (4 Positions)
- Height range of 47cm to 101cm
- Tilting angle from horizontal to +90°
- Lowers to wheelchair height for ease of patient transfer
- Tilt angle inclinometer fitted as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Supplied with patient hand grips and set of harnesses
- Worktable supplied as standard
- 70cm wide upholstery as standard

Optional Accessories

- 6057** Battery Back-Up
- 6058** Lateral Supports (Pair)
- 6059** Adjustable Pommel Support*
- 6007** Head Support Cushion
- 6062** Standard Harness Set
- 6063** Dual Leg Harness Set*
- 6064** Pair Dual Knee Harnesses*
- 6065** Standard Harness (Wide)
- 6066** Standard Harness (Narrow)

* (will fit ST7641DL only).



Central locking castors



Angle inclinometer



Model ST7641DL Divided leg model with optional pommel and lateral supports. Colour light grey