



\*When fitted with 9202 Easy Clean Upholstery

## Therapy 2 Section Basic Head Couch



Designed for head and neck mobilisation, as well as manual traction techniques, these strong and rigid tables feature a short head section, with the option of adjustable height armrests on the "Plus" model, providing additional patient comfort during prone treatments. The head section has a breathing hole and plug fitted as standard.

### Models

- ST2558 Hydraulic with manual backrest
- ST2568 Electric with manual backrest



Model ST2558 Hydraulic model with manual backrest. Colour mandarine

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- The power assisted head section incorporates anti-trap design gas strut
- Head section angle adjustment from -25° to +40°
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery

### Optional Accessories

- 6001 Perimeter Foot Switch
- 6200 Single Foot Switch
- 6239 Dual Foot Switch Operation
- 6006 Upper Limb Support
- 6002 Side Support Rails
- 6002P Side Support Rail Cushions
- 6005 Paper Roll Holder
- 6007 Removable Head/Knee Support
- 6053 Protective Base Shroud
- 6010 Couch Cover (Breathing Hole)
- 6011 Couch Cover (Without Breathing Hole)

See page 48-55 for additional information on our range of optional accessories and colour options



Negative head section to -25°



Low patient transfer height