



Therapy 2 Section Plus Head Couch



The plus head design is ideal for neck mobilisation as well as manual traction techniques. These couches feature a short head section with height adjustable armrests for prone treatments, promoting relaxed shoulders and a more comfortable patient. The couch has an excellent height range of 45cm to 98cm and is fitted with a breathing hole and plug as standard.

Models

- ST2559** Hydraulic with manual head section
- ST2569** Electric with manual head section



Model ST2569 Plus head model with manual head section and optional perimeter foot switch. Colour tileul

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- The power assisted head section incorporates anti-trap design gas strut
- Head section angle adjustment from -25° to +40°
- Adjustable height armrests as standard
- Choice of manual or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery

Optional Accessories

- 6001** Perimeter Foot Switch
- 6200** Single Foot Switch
- 6239** Dual Foot Switch Operation
- 6006** Upper Limb Support
- 6002** Side Support Rails
- 6002P** Side Support Rail Cushions
- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6012** Couch Cover (Breathing Hole) for Plus Head Models (4 Piece)

See page 48-55 for additional information on our range of optional accessories and colour options



Positive head section to +40°



Low patient transfer height