



*When fitted with 9201 Easy Clean Upholstery

Therapy 2 Section Standard Head Couch



A very strong and rigid 2 section couch with a wide range of optional accessories available and an extensive choice of colours. This couch is suitable for supporting very heavy patients for examinations, minor procedures and physiotherapy treatments. This 2 section couch is available as a hydraulic or electric model, both models feature a manual backrest that can be adjusted from -25° to +75° and has an excellent height range from 45cm to 98cm.

Models

- ST2557** Hydraulic with manual backrest
- ST2567** Electric with manual backrest



Model ST2567 Electric model with manual backrest and optional foot switch. Colour light grey

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- The power assisted backrest section incorporates anti-trap design gas strut
- Backrest angle adjustment from -25° to +75°
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard upholstery width of 65cm

Optional Accessories

- 6001** Perimeter Foot Switch
- 6200** Single Foot Switch
- 6239** Dual Foot Switch Operation
- 6006** Upper Limb Support
- 6002** Side Support Rails
- 6002P** Side Support Rail Cushions
- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6010** Couch Cover (Breathing Hole)
- 6011** Couch Cover (Without Breathing Hole)

See page 48-55 for additional information on our range of optional accessories and colour options



Negative head section to -25°



Low patient transfer height