



*When fitted with 9302 Easy Clean Upholstery

Therapy 3 Section Basic Head Couch



Designed with a short head section for head and neck mobilisation techniques, these versatile tables are suitable for a wide range of physiotherapy and rehabilitation procedures. The couch has an excellent height range of 45cm to 98cm and is fitted with a breathing hole and plug as standard.

Models

- ST3558** Hydraulic with manual backrest and head section
- ST3568** Electric with manual backrest and head section



Model ST3558 in seated position.
Colour light grey

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Power assisted head section is adjustable from -25° to +40°
- Breathing hole and plug fitted as standard
- Power assisted backrest is adjustable from 0° to +80°
- Head section and backrest incorporate anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery

Optional Accessories

- 6001** Perimeter Foot Switch
- 6200** Single Foot Switch
- 6239** Dual Foot Switch Operation
- 6006** Upper Limb Support
- 6003** Side Support Rails
- 6003P** Side Support Rail Cushions
- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6010** Couch Cover (Breathing Hole)
- 6011** Couch Cover (Without Breathing Hole)

See page 48-55 for additional information on our range of optional accessories and colour options



Negative head/foot section to -25°