



Therapy 3 Section Plus Head Couch



Designed with a short head section for head and neck mobilisation techniques, these versatile tables are suitable for a wide range of physiotherapy and rehabilitation procedures. The "Plus" model features adjustable height armrests which promote a relaxed shoulder girdle and a more comfortable patient during prone treatments. The couch has an excellent height range of 45cm to 98cm and is fitted with a breathing hole and plug as standard.

Models

- ST3559** Hydraulic with manual backrest
- ST3569** Electric with manual backrest

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Power assisted head section is adjustable from -25° to +40°
- Adjustable height armrests for prone treatments
- Power assisted backrest is adjustable from 0° to +80°
- Head section and backrest incorporates anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery

Optional Accessories

- 6001** Perimeter Foot Switch
- 6200** Single Foot Switch
- 6239** Dual Foot Switch Operation
- 6006** Upper Limb Support
- 6002** Side Support Rails
- 6002P** Side Support Rail Cushions
- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6012** Couch Cover (Breathing Hole) for Plus Head Models (4 Piece)

See page 48-55 for additional information on our range of optional accessories and colour options



Model ST3569 in seated position with optional perimeter foot switch.
Colour cobalt blue



Model ST3569 at mid-height with optional perimeter foot switch.