



Therapy 3 Section Standard Head Couch



A very strong and rigid 3 section couch with a traditionally long head section to suit a wide range of physiotherapy and general examination procedures. The couch can be adjusted to provide support for the patient in prone, supine and seated positions. The couch has an excellent height range of 45cm to 98cm and is fitted with a breathing hole and plug as standard.

Models

- ST3557** Hydraulic with manual backrest and head section
- ST3567** Electric with manual backrest and head section



Hydraulic foot operation



Model ST3557 Colour sky blue



Negative head/foot section to -25°



Low patient transfer height

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Power assisted head section is adjustable from -25° to +75°
- Power assisted backrest adjustable from horizontal to 80°
- Head section and backrest incorporate anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery

Optional Accessories

- 6001** Perimeter Foot Switch
- 6200** Single Foot Switch
- 6239** Dual Foot Switch Operation
- 6006** Upper Limb Support
- 6003** Side Support Rails
- 6003P** Side Support Rail Cushions
- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6010** Couch Cover (Breathing Hole)
- 6011** Couch Cover (Without Breathing Hole)

See page 48-55 for additional information on our range of optional accessories and colour options