



ST4552
ST4562
ST4542



ST4552W
ST4562W
ST4542W



Therapy Bobath Couch



By tradition, Bobath treatment techniques require a suitably wide treatment surface with a safe and capable lifting capacity for both patient and therapist. The Bobath Therapy range of couches has more than achieved this with a choice of 2 sizes and higher than standard safe working load of 250Kg. See deluxe model for 325Kg safe working load.

Models

ST4552	Hydraulic with manual backrest 105cm
ST4562	Electric with manual backrest 105cm
ST4542	Electric with electric backrest 105cm
ST4552W	Hydraulic with manual backrest 125cm
ST4562W	Electric with manual backrest 125cm
ST4542W	Electric with electric backrest 125cm



Model ST4562 Shown with optional perimeter foot switch. Colour dark blue

Standard Features

- Safe working load and lifting capacity of 250Kg (550lbs)
- The power assisted backrest section incorporates anti-trap design gas strut
- Backrest angle adjustment from horizontal to +75°
- Choice of manual or electric height elevation
- Excellent height range from 46cm to 100cm
- Retractable wheel design with adjustable foot for stability on uneven floors and surfaces
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Choice of widths available 105cm or 125cm

Optional Accessories

- 6004 Breathing Hole & Plug
- 6057 Battery Back-Up
- 6207 Double Foot Switch for Models ST4542/W
- 6001XL Perimeter Foot Switch
- 6074 Side Support Rails
- 6074P Side Support Rail Cushions

See page 48-55 for additional information on our range of optional accessories and colour options

Model ST4562 Shown with optional perimeter foot switch in low patient transfer height.

