



Therapy Deluxe Drainage Couch



The deluxe drainage couches have the same head section options as the standard 2 and 3 section therapy couches and also provides the option of having the backrest/foot section supplied in a split leg model for additional mobilisation techniques. The electric drainage facility can be operated on either side, via a simple push button hand switch mounted near the top frame. This facility provides angular adjustment to +25° at the hinge point of the of the mid and foot sections. The perimeter foot switch system is standard on these models.

Models

ST3547	Standard head design
ST3547S	Standard head design with split leg
ST3548	Basic head design
ST3548S	Basic head design with split leg
ST3549	Plus head design
ST3549S	Plus head design with split leg

Model ST3549S with plus head design and split leg. Colour dark blue



Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric drainage facility to +25° between mid section and backrest
- Electric height operation via integral perimeter foot switch
- Excellent height range from 46cm to 99cm
- Adjustable angle backrest from horizontal to +90°
- Adjustable head section with positive and negative angle
- Power assisted backrest and head sections incorporate anti-trap design
- Choice of standard, basic or plus head design with adjustable armrests
- Choice of standard or split leg foot section
- Breathing hole and plug fitted as standard
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Standard width 65cm upholstery

Optional Accessories

- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6100** Matching Operators Chair



Couch in horizontal position