

Width
65cm



130cm

Small

Width
65cm



152cm


Medium

Width
65cm



186cm

Large



Height range 41cm to 94cm

Therapy Hygiene Tables



For environments such as special schools, where many of the children require regular hygiene attention, these variable height tables provide not only a low transfer height but quick and simple adjustment to the correct height for changing procedures, reducing back and postural problems for support staff. Standard width 65cm. Wider models also available. All large models are 186cm long.

186cm Models

- ST1551A** Hydraulic with retractable wheels
- ST1561A** Electric with retractable wheels
- ST1651A** Hydraulic with central locking castors
- ST1661A** Electric with central locking castors



Model ST1551A Shown with optional side support rails.



Model ST1551A Large hydraulic model, shown with retractable wheels and optional side support rails and cushions. Colour lotus green

152cm Models

- ST1551B** Hydraulic with retractable wheels
- ST1561B** Electric with retractable wheels
- ST1651B** Hydraulic with central locking castors
- ST1661B** Electric with central locking castors

130cm Models

- ST1551C** Hydraulic with retractable wheels
- ST1561C** Electric with retractable wheels
- ST1651C** Hydraulic with central locking castors
- ST1661C** Electric with central locking castors



Model ST1551C low transfer height. Colour zest



Model ST1551C Small hydraulic model with retractable wheel design. Colour zest



Central locking castors