

SOL

Rotating podiatry chair equipped with a hydraulic pump that controls height. Backrest and legrests are independently adjustable by gas lift. White PU seamless upholstery of high quality that eases cleaning.

Sillón de podología rotatorio. Altura regulable mediante bomba hidráulica. Respaldo y reposapiés regulables independientemente mediante pistón de gas. Tapizado blanco en PU sin costuras de alta calidad y fácil de limpiar.

Fauteuil de podologie rotatif. Hauteur réglable par pompe hydraulique. Dossier et jambières réglables de façon indépendante par vérin à gaz. Revêtement en PU blanc sans coutures de haute qualité et facile à nettoyer.

Drehbarer Fußpflegestuhl. Die Höhe ist mittels einer Hydraulikpumpe verstellbar, die Rückenlehne und die zwei unabhängigen Fußstützen sind mittels Gasfeder regulierbar. Hochwertige weiße PU-Polsterung ohne Nähte, leicht abwaschbar.

Poltrona podologica girevole regolabile in altezza con una pompa idraulica. Schienale e gambali regolabili indipendentemente con un pistone a gas. Rivestimento in PU bianco di alta qualità senza cuciture e facile da pulire.

Podiatry Hydraulic Chair Ref. 2231A.A26

Size without armrest: 182x56x56/67 cm

Size with armrest: 182x74,5x56/67 cm

Weight: 37 Kg

Upholstery: PU

- 1** Angle 45° ~ 85°
- 2** Rotation 360°
- 3** Angle 0° ~ 85°
- 4** Extension ±15
- 5** Up and Down 56/67 cm

OPTIONS

Bed cover

Paper roll holder

6 *Debris tray* Ref. 30030

PACKING SIZE cm

140x70x70

N.W. kg

37

G.W. kg

46

MAX WEIGHT

150 kg

6



1

3

4

2

5



SOL

