

# FINEX

## Digital paraffin heater for hand and foot care with timer. Capacity: 5 l

Fundidor de parafina digital para tratamiento de manos y pies. Equipado con temporizador. Capacidad: 5 l

Chauffe-paraffine digital pour le soin des mains et pieds. Il est équipé d'un temporisateur. Capacité: 5 l.

*Digitale Paraffinerhitzer für Hände und Füße mit Timer ausgestattet. Füllmenge: 5 L.*

*Scaldaparaffina digitale per trattamenti di mani e piedi. Dotato di un temporizzatore. Capacità: 5 l.*

## Digital Paraffin Heater Ref. WKE004

Size: 44,5x27x19 cm  
Weight: 3Kg

Input Voltage: 220~240V  
Rated power: 220W  
Temperature range: 46-60°C  
Complete melting time: 30 minutes  
Capacity: 5l.

### PACKING SIZE cm

48x29x22

N.W. kg

3

G.W. kg

3,7

